

Instructions for use

MASTER lase EXPERT lase



English

C€ 0124

Order number for these instructions for use: 1.010.0274



© 2012 - 2013 Kaltenbach & Voigt GmbH

These instructions for use are addressed to doctors, dentists and medical assistants who have been trained in the intended use of this medical device in accordance with the statutory directives.

The duplication of these instructions for use, also extracts thereof, is only permissible with the approval of Kaltenbach & Voigt GmbH.

All rights shall be reserved for the case of patent registration or registration of utility model.

All information contained in these instructions for use have been checked to the best of our knowledge and judgement. However, neither author nor publisher can be held liable for any damage arising in association with the application of these instructions for use.

Names of companies, products or persons listed in this document as well as the characters and/or data included are ficticious and in no circumstances represent any real, existing individual person or existing company, product or outcome unless manufacturer related.

Kaltenbach & Voigt GmbH is in possession of numerous patents, patent applications, trademarks, copyrights or other froms of intellectual property applicable to this document. In as far as this had not been explicitly approved in a written licensing agreement by Kaltenbach & Voigt GmbH, the user or owner shall not have any rights to these patents, patent applications, trademarks or other intellectual property with the provision of this document.

Other products, brands and company names listed in this document could be the trademarks of their respective owners. The adherence to all applicable copyright laws is the responsibility of the user of this document.



Trademarks

Registered trademarks, trade names and names are used in these instructions for use. The relevant protective provisions are applicable even if they are not designated as such.

Manufacturer



Kaltenbach & Voigt GmbH

Bismarckring 39

D-88400 Biberach/Riß

Service Support

KaVo Dental GmbH

Bahnhofstraße 20

88447 Warthausen

Phone: +49 7351 56-2800

Mail: Service.Laser@kavo.com



Table of contents

	Table of contents	1
	Use of these instructions for Use	5
lr	ntroduction6	;
	Intended medical use - application field	6
	Indications, indication areas	7
	Contraindications	7
	Requirements for laser application	7
	Education and declaration of consent of patient	7
S	afety8	}
	Medical Device Law (MPG) / EU – Medical products - guideline	8
	Possible complications, risks and hazards	9
	Installation and initial commissioning	. 10
	Laser Safety	. 11
	Application of laser radiation	. 13
	Laser safety office	. 13
	Eye protection	. 14
	Safety of the patients, user and third parties	. 15
	Training and responsibility of the user	. 15
	Training and 75B due dilligence of the operating personnel	. 16
	Responsibility of Kaltenbach & Voigt GmbH	. 16
	Warning and safety instructions concerning electromagnetic compatibility accordance with EMC - Norm DIN EN 60601-1-2:2007	y in . 17
	Guidelines and manufacturer's declaration - electromagnetic emission	ւ:17
	Guidelines and manufacturer's statement - Electromagnetic immunity:	18
	Safety measures and safety systems	. 22
	Electrical installation	. 22
	Disconnection from the supply network	. 22
	Laser Stop	. 22
	Laser safety chain	. 23



Laser audible warning	23
Safety during transportation	23
Operating position	23
Exchanging parts	23
Combination with accessories	24
Sterilisation of accessories	24
Warning and instruction signs	25
Description	28
Technical specifications	28
Scope of delivery and accessories	30
Setup	36
Terminals	37
Electrical terminals	37
Transmission optic fibres connection	37
Operating elements of the wireless foot switch	38
Power supply of the wireless foot switch	39
Power supply of the battery	40
Control panel	41
Therapy mode	43
Setting parameters	43
Selection of indication	45
Pilot	48
Protocol	48
Menu	49
Operation	51
Observe safety	51
Preparation for commissioning	52
Mounting the charger adaptor	53
Dismounting the charger adaptor	
Assemble wireless foot switch and log into the device	55



С	ommissioning	. 57
	Commissioning of the transmission fibre optic conductor	. 59
	Mounting the sapphire tip handpiece and commissioning	. 62
	Mounting the bare fibre handpiece and commissioning	. 66
	Blacken the sapphire tip / bare fibre tip	. 70
	Dismounting the sapphire tip handpiece	.72
	Dismantling the bare fibre handpiece	. 73
Т	herapy settings	. 74
	Select the indication	.74
	Setting the performance	. 75
	Setting the performance with the wireless foot switch	. 76
	Setting the pulse form	. 77
	Setting the pulse duration	. 78
	Setting the pulse pause	. 79
	Setting the duration of application	. 80
	Compilation and use of list of Favorites	. 81
	Saving the default settings of the indication parameters	. 81
	Saving individual indication parameter settings	
	Calling up indication parameters from the list of Favorites	
	Deleting indication parameters from the list of Favorites	
	Changing the display mode of the protocol button	
	Resetting the protocol or timer	
N/	lenu settings	
IV	Setting the volume for laser warning signal tones and input feedback.	
	Setting the switch-off time	
	<u> </u>	
	Setting transitional period Ready → Standby	
	Selecting option "Pilot during treatment"	
	Selecting operation	
	Setting the language	
	Setting the date and time	. 95

MASTER lase / EXPERT lase



	Changing the PIN	96
	Setting the frequency of the PIN - request	98
	Reading the logbook	98
	Triggering the therapy beam	100
	Triggering the therapy beam with the wired foot switch	101
	Triggering the therapy beam with the wireless foot switch	102
	Switching off the device	103
	Transportation of the device	104
Ma	aintenance, cleaning, safety check	105
	Maintenance by the operator	105
	Classing	100
	Cleaning	106
	Safety checks	
,	-	108
M	Safety checks	108 109



Use of these instructions for Use



Warning means that if instructions are ignored or not followed correctly, this could lead to the endangering of patients, and/or operator and/or third parties.



Caution means that if instructions are ignored or not followed correctly, this could lead to damage of the device.



Instructions are given by additional information.

A cross reference with item number is provided in brackets: e. g. (4, fig. 3)

meaning: the component mentioned is presented in fig. 3, item 4.

Operating instructions are marked with bold numbers.

e.g. 1 Switch the device on.

2 Switch peripherals on.

A single operating instruction is also assigned a number.

If a setting is made by directly touching a display, the symbol is used in combination with the display text: e. g. *Pilot*-settings on the pilot light are selected by touching the display text "Pilot".

Display indications are presented in italics: e. g. WATT



Images and on-screen displays in the present instructions for use may be slightly different in appearance from the actual product as long as this does not have any adverse effects on the technical features and safety of the product.



Introduction

Intended medical use - application field

The device is a portable diode laser for multidisciplinary application in dental surgeries and practices.

The power transmitted to the tissue by the device amounts to up to 10 Watt.



The device may only be used in accordance with its intended use.

The manufacturer shall assume no liability in case of application not conformant to intended use.



Use or operation that does not comply with the description in these instructions for use could lead to dangerous exposure to radiation.

The manufacture shall assume no liability for interventions in the product.



Alterations to the device are prohibited.



Indications, indication areas

The indication area of the MASTER/EXPERT lase is subject to continual extension due to clinical research. State-of-the-art findings can be found in the current specialised literature and current publications.

The MASTER/EXPERT lase can be employed for dental applications within the following indication areas:

- Dental-surgical tissue preparations such as resection, extirpation, ablation, vaporisation, coagulation, haemostasis, incision, excision;
- Germ reduction:
- Desensitisation, bio-stimulation;
- Activation of bleaching material.

Dental-surgical indications are primarilly given in operations, in which a minimum invasive intervention, good haemostasis, extended coagulation volume or cutting and removal of tissue with the least possible blood loss is in the foreground.

Indications for selection can be found in table: Selection of indications, page 46.

Contraindications

Those patients who suffer from photodermatosis and who are photosensitive (photoallergies) may not be treated.

Requirements for laser application

The fundamental prerequisite for successful, safe application of the laser for medical reasons is a comprehensive specialised knowledge of laser knowledge.

Current addresses and notes on laser literature can be requested from Kaltenbach & Voigt GmbH.

Experimental work helps quickly and reliably to understand the connections between laser radiation and tissue samples under different treatments modes and parameters, types of tissue and different accessories.

Education and declaration of consent of patient

The patient must be educated about the course of treatment. The customary declaration of consent must be obtained from the patient.



Safety

Prior to the startup of the MASTER/EXPERT lase, every user must have read the complete instructions for use and understand them.

The adherence to the safety guidelines when handling the laser is the responsibility of the laser safety officer.

All safety aspects defined in the section "Safety" of these instructions for use must be adhered to.

The instructions for use must always be available in the vicinity of the device.

The instructions for use must be passed on to all subsequent users.

Medical Device Law (MPG) / EU - Medical products - guideline

The laser device is an active medical device in accordance with Class IIb pursuant to Directive 93/42/ECC concerning medical devices.

In accordance with the Medical Product Law, medical devices may only be set up, used and operated in accordance with their intended purpose, according to the guidelines of this law and the pursuant approved legislation, the acknowledged rule of technology and the occupational safety and accident prevention regulations. They may not be used or operated if they exhibit defects or could endanger patients, staff or other persons. They may only be used and operated by persons who can ensure proper handling based on their training or knowledge and practical experience.

Class IIb medical devices may only be used by persons who have been familiarised with the device and trained in appropriate handling in consideration of the instructions for use.

Only such persons who can ensure proper handling based on their training or knowledge and practical experience in handling the device may be trained.

In addition, the user must keep a medical device log.



Possible complications, risks and hazards

Independently of laser, all known operative and clinical risks and complications must be taken into account which are due to the operative procedure, the condition of the patient and the clinical pattern.

There a laser-specific risks and hazards for the patient or surgical personnel due to:

- Thermal damage of such critical body structures that are in the immediate vicinity of the treatment area (e.g. teeth. bone, periodontium).
- Thermal damage to critical body structures due to unintentional activation and the associated remote effect of the laser beam
- Damage to sensitive body parts such as the eye of skin of the operating personnel or patients due to unintentional, undirected activation of the laser beam in the laser area.

⚠ WARNING

Potential interaction of the laser light with optical diagnostic and monitoring systems.

In case of intended use of high intensity light, the device transmits such wavelengths that are defined in the laser specification label of the individual device. Attention must be paid both to the wavelength of the therapy laser and of the pilot laser. Irradiation of this light (direct or due to scattering) in the detectors of optical diagnosis and monitoring systems can result in measuring errors or faulty interpretation of measurements by this device. This can lead to potentially dangerous situations for the patient.

Avoid irradiation of the light, both of the therapy laser and the pilot laser (direct or due to scattering) in the detectors of optical diagnosis and monitoring systems.



Installation and initial commissioning

The installation and initial commissioning may only be performed by Kaltenbach & Voigt GmbH staff or by persons authorised by Kaltenbach & Voigt GmbH. The devices may not be used until:

- acceptance has been obtained from the purchasing body
- a laser protection officer has been appointed by the purchasing body or the operator*
- the device has been registered with the German Federation of Institutions for Statutory Accident Insurance and Prevention and the authorities responsible for occupational safety*
- All safety measures have been carried out
- the device is mastered, even in cases of emergency, due to training of the personnel
- the device has been handed over ready for operation by Kaltenbach & Voigt GmbH and the personnel has been trained

The inventory list and instructions for use must be kept in an accessible location.

Accidents and personal injuries must be reported to Kaltenbach & Voigt GmbH or the responsible authority without delay.

^{*} Applicable if required by national industrial standards.



Laser Safety

All valid guidelines concerning radiation protection for laser devices must be adhered to when operating the device. In case of ambiguity of the guidelines, the laser safety officer must be consulted.



Irreversible damage could be caused.

The retina and the skin may not be exposed to direct or reflected laser radiation (e.g. due to shining materials) under any circumstances. The laser safety glasses provide only short-time protection from direct laser light.

The laser area is defined as the area in which the values for the maximum permissible radiation dose (MRD) for eyes (IEC/EN 60825-1:2007) could be achieved or exceeded.

The entire operating room is regarded as laser area if hand applicators are used for open work.

An unintentional distraction of the laser beam must be taken into consideration.

Adequate laser safety goggles must be worn inside the laser area.

The laser area must be marked with door signs and warning lights, must be kept as small as possible and secured against unauthorised entrance.

The number of persons within the laser area must be as low as possible.

Reflecting, shining materials must be removed form the laser area or packed in flame-retardant cloths. Combustible materials must be removed.



The application of laser radiation can lead to smoke formation. This can contain viable particles of tissue.

Provide sufficient smoke extraction in the laser area.

Substances, which could contribute to the generation of gases, dust and mist due to the effect of laser radiation, as well as inflammable substances, must be removed from the laser area or be protected by appropriate measures.

Only devices that are suitable for laser treatment may be used during laser treatment. Special laser devices inhibit the reflection of the laser beam by their shape and surface characteristics.

Appropriate protective equipment must be provided in adequate quantity and flawless condition.



If endoscopes are used without a video camera, appropriate laser protection filters must be used for the endoscope. The user need not wear safety goggles during laser endoscopy with laser filter. All other persons present must wear safety goggles.

The laser beam must not be used in body regions with oxygenenriched air or where combustible gases or vapours exist or could arise. It is essential that protective measures are taken against fire and explosion hazards.



The device generates a highly energetic laser ray.

The device may not be operated in potentially explosive and/or burns inductive atmospheres.

An explosive or inflammable atmosphere is deemed to exist e.g. on incidence of vapours arising from anaethesia, cleaning agents or disinfectants, or due to oxygen enrichment as a result ventilation using oxygen.

All persons working in the laser area must receive annual training on laser safety and the operation of the devices. Participation in training must be confirmed in writing.



Application of laser radiation

The device is a laser device of Class 4 according to IEC 60825-1:2007 / EN 60825-1:2007.



High-energy laser beam.

Inappropriate use can lead to personal injury and/or damage to devices.

All prescribed safety measures must be accomplished and the corresponding safety aids used and kept in faultless condition.

Laser safety office*

According to the Laser Safety Directive BGV B2 (Germany), a laser safety officer must be appointed in writing by the operator or owner of the equipment. The laser safety officer is responsible for:

- Performance and verification of the safety measures and the protective equipment
- Familiarisation of the users with the safety measures and operation of the device
- Identification of the laser area
- Checking the warning lights at the entrance to the laser area
- Safe keeping of the access code (PIN)
- Safe storage of the device
- Correct connection of the device following a change of location

24.06.2013 GA_MASTERlase+EXPERTlase_1.010.0274_03_en.docx

^{*} Applicable if required by national industrial standards.



Eye protection

The eye is especially endangered due to the high energy density. The eye can be damaged even by weak laser radiation.

Laser area according to IEC EN 60825 (version 03-2007; section 3.61).

The safe distance from the eye to the laser beam egress point is defined as "Nominal Ocular Hazard Distance" (NOHD). Laser safety goggles need not be worn outside this area if the fibre optic conductor is not operated in combination with the focussing optics. The beam divergence (a) is 10° (full angle).

The safety distance NOHD for the MASTER/EXPERT lase is defined as 3.2 m to the distal end of the fibre optic conductor or handpiece used.



Irreversible eye damage could be caused.

Each person within the laser area must wear safety goggles.

Laser safety goggles for the laser system must fulfil the following conditions at least:

Laser type	Wavelength of the filtered beam	Protection class EN 207:2002	Manufacturer abbreviation
D	980 nm	LB3	in accordance with EN 207:2009

Laser safety goggles do not provide absolute protection, therefore the exposure time should be brief. The safety goggles may not exhibit any damage.



The recommended protection class is only applicable for the use of accessories without imaging optics, e.g. fibre optic conductor with blank fibre tip. An effective protection for accessories with imaging optics, e.g. focussing handpieces, is only guaranteed when they are used as intended. Furthermore, the relevant national industrial standards in their valid version must be adhered to, e.g. IEC 60825 (Safety of Laser Devices)



Safety of the patients, user and third parties

Numerous safety devices are incorporated into the device.

All functions are continuously monitored automatically.

The therapy laser can only be activated if all functions are faultless.

Training and responsibility of the user

As does every highly effective medical device, a laser device demands special expertise and caution in handling and use. The laser device may only be used by persons who have been familiarised with the device and trained in appropriate handling in consideration of the instructions for use. and are familiar its therapeutic effects and potential hazards.

Training courses are conducted by specialists authorised by Kaltenbach & Voigt GmbH as well as by reference physicians with clinical expertise.



Irreversible damage to health could be caused.
Read and observe the above mentioned information.
The safe and successful performance of the therapy requires participation in a training course prior the initial use.

Prior to every treatment, the user must have verified the proper functioning of the device in order to exclude any danger to patients, third parties or himself. If the device is used with accessories, the instructions for use of the accessories must have been read completely and understood.



Training and due dilligence of the operating personnel

Untrained or unqualified operating personnel may not operate the device in any circumstances.

Training courses are conducted by specialists authorised by Kaltenbach & Voigt GmbH.

Operating and cleaning personnel must exercise enhanced care when handling the device. More information in the chapter on cleaning.

Responsibility of Kaltenbach & Voigt GmbH

In relation to the safety, reliability and performance of the device, Kaltenbach & Voigt GmbH shall only be held responsible if:

- The installation, adjustment and initial commissioning are only performed by Kaltenbach & Voigt GmbH staff or by persons authorised by Kaltenbach & Voigt GmbH
- The electrical system in the relevant room meets the requirements of the national standards
- The unit is operated according to the instructions for use



Warning and safety instructions concerning electromagnetic compatibility in accordance with EMC - Norm DIN EN 60601-1-2:2007

Medical electrical devices are subject to special measures regarding electromagnetic compatibility and must be operated, installed and commissioned in accordance with the EMC instructions contained in the KaVo instructions for use.

Portable and mobile high-frequency communications devices may interfere with electrical medical devices.

The EMC requirements according to EN 60601-1-2:2007 are applicable in combination with the power cable connected to the device, type H05VV (length max. 2.5 m)



Using other accessories or cables as those defined - with the exception of internal original spare parts - can lead to increased emission or reduced interference immunity in the device.

Guidelines and manufacturer's declaration - electromagnetic emission:

The device is for use in an electromagnetic environment like the one cited below. The customer or user of the device should ensure that it is used in the correct environment.

Measurement of emissions	Conformance	Electromagnetic environment - Guidelines
HF emissions according to CISPR 11	Group 1	The device uses HF energy for its internal functions exclusively. Its HF transmission is therefore very low, and it is improbable that neighbouring electronic devices will be disturbed.
HF emissions according to CISPR 11	Class B	The device is for use in all facilities including residential areas and facilities that are directly connected to a public power supply that also supplies
Harmonics in accordance with IEC 61000-3-2	Class A	residential buildings.
Voltage fluctuation / flicker in accordance with IEC 61000-3-3	complies	



⚠ CAUTION

The device may not be used or stacked directly next to other devices.

If operation of the device stacked right next to other devices is indispensable, the device must be monitored to verify the proper operation in this arrangement.

Guidelines and manufacturer's statement - Electromagnetic immunity:

The device is for use in an electromagnetic environment like the one cited below. The customer or user of the device should ensure that it is used in the correct environment.

Immunity tests	IEC 60601 test levels	Compliance level	Electromagnetic environment - Guidelines	
Electrostatic discharge (ESD)	± 6 kV contact discharge	± 6 kV contact discharge	Floors should be made of wood or concrete or be fitted with ceramic tiles. If	
according to IEC 61000-4-2	± 8 kV atmospheric discharge ± 8 kV atmospheric discharge		the floor is made of synthetic material, the relative humidity must be at least 30%.	
Fast transient electrical	± 2 kV for power lines	± 2 kV for power lines	The quality of the supply voltage should correspond to that of a typical business or hospital environment.	
disturbances/ bursts according to IEC 61000-4-4	± 1 kV for input and output lines	± 1 kV for input and output lines		
Surges according to	± 1 kV push-pull voltage	± 1 kV push-pull voltage	The quality of the supply voltage should correspond to	
IEC 61000-4-5	± 2 kV common mode voltage	± 2 kV common mode voltage	that of a typical business or hospital environment.	



Immunity tests	IEC 60601 test levels	Compliance level	Electromagnetic environment - Guidelines	
	< 5 % U _{T for ½ period} (> 95 % interruption)	< 5 % U _{T for ½ period} (> 95 % interruption)	The quality of the supply voltage should correspond to that of a typical	
Voltage interruptions,	40 % U _{T for 5 periods} (60 % interruption)	40 % U _{T for 5 periods} (60 % interruption)	business or hospital environment.	
short-term interruptions, and fluctuations of the supply voltage	70 % U _{T for 25 periods} (30 % interruption)	70 % U _{T for 25 periods} (30 % interruption)	When the user of the device needs continued operation even when the power supply is interrupted, it is recommended to supply the device from an uninterrupted power supply or a battery.	
according to IEC 61000-4-11	> 5 % U _{T for 5 s} (> 95 % interruption)	> 5 % U _{T for 5 s} (> 95 % interruption)		
Magnetic field at a supply frequency (50/60 Hz) according to IEC 61000-4-8 mains freque should corres the typical value a business at hospital		Magnetic fields at the mains frequency should correspond to the typical values in a business and hospital environment.		
NOTE: U _T is the alternating mains voltage before the test level is used				



Immunity tests	IEC 60601 test levels	Compliance level	Electromagnetic environment - Guidelines		
including cables than th	Handheld and mobile wireless devices should not be used at a shorter distance from the including cables than the recommended safe clearance calculated with the appropriate equation for the emission frequency.				
			Recommended safe distance:		
Conducted HF disturbances according to IEC 61000-4-6	10 V _{eff} 150 kHz to 80MHz in the ISM bands ^a	1 V _{eff}	$d=1.2 \sqrt{P}$		
Wireless HF interference according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	d = 1.2 \sqrt{P} 80 MHz to 800 MHz		
			d = 2.3. \sqrt{P} 800 MHz to 2.5 GHz		

With P as the rated power of the transmitter in Watts (W) according to the transmitter manufacturer, and d as the recommended safe distance in meters (m).

The field strength of stationary wireless radio transmitters as measured locally:

b should be lower than the conformance level of all frequencies.

^c interference is possible in the vicinity of devices that bear the following symbol.



NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE 2:These guidelines may not be applicable in all situations. The spread of electromagnetic waves is absorbed and reflected by buildings, objects and people.

^a The ISM bands between 150 kHZ and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.

^b The field strength of stationary transmitters such as base stations of mobile telephones and land radio devices, amateur radio stations, AM and FM, radio and television broadcasters cannot be theoretically predetermined accurately. A site study should be considered to determine the electromagnetic environment in terms of stationary transmitters. The device should be observed to verify the functioning as intended. Should unusual performance features be observed, additional measures may be required, such as e.g. a different alignment or another location for the device.

 $^{\rm c}$ The field strength is less than 3 V/m over the frequency range of 150 kHz to 80 MHz.



Recommended safe distance between portable and mobile HF telecommunications equipment and the device

The device is intended for use in an electromagnetic environment in which HF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic disturbances by maintaining the minimum distance between portable and mobile HF telecommunications devices (transmitters) and the device - depending on the output of the communication device as indicated below.

Nominal power of the transmitter (W)	Safe distance according to the transmission frequency (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=$ 1.2 \sqrt{P}	d = 1.2 \sqrt{P}	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose maximum rated power is not in the above table, the distance can be calculated using the equation for the respective gap, where P is the maximum rated power of the transmitter in Watts (W) according to the manufacturer's information.

NOTE 1:

To calculate the recommended safe distance from transmitters with a frequency range of 80 MHz to 2.5 GHz, an additional factor of 10/3 was used to reduce the probability that a mobile/portable communication device unintentionally brought into the patient area would cause malfunction.

NOTE 2:

These guidelines may not be applicable in all situations. The spread of electromagnetic waves is absorbed and reflected by buildings, objects and people. A site study should be considered to determine the electromagnetic environment in terms of mobile transmitters. The device should be observed to verify the functioning as intended. Should unusual performance features be observed, additional measures may be required, such as e.g. a different alignment or another location for the device.



Safety measures and safety systems

Electrical installation

The device may only be used in rooms designed for medical purposes in which installations are carried out in accordance with the relevant national installation standards. The device may only be connected to a separate secured socket. The device is a protection class 2 device according to IEC/EN 60601-1 (revision 1996 & 2005).

For proper operation in relation to the electromagnetic compatibility, see "Warning and safety instructions concerning electromagnetic compatibility in accordance with EMC - Norm DIN EN 60601-1-2:2007".

The distance between the edges of the housing of the device and the wall must be at least 10 cm to avoid impairing ventilation and exhaust. Installation in closed cabinets in not permissible.

Disconnection from the supply network

Pull the power plug to disconnect the device from the power supply. Ensure unobstructed access to the power plug when installing the device.

Laser Stop

A LASER STOP button is provided to enable quick switching off by the operator during operation.

Figure 1. Laser STOP button






Laser safety chain

A chain of safety devices is incorporated into the device. This laser safety chain reliably prevents unintentional triggering. The laser safety chain consists of hardware switches arranged in sequence.

The laser cannot be triggered if a switch is open, i.e. a safety condition is not fulfilled. The final link of the laser safety chain is the shutter, an electronic switch which short-circuits the laser diode, thus disconnecting the power supply.

Figure 2. Laser safety chain



Laser audible warning

During the laser activation, an acoustic signal is emitted to warn of laser radiation.

Safety during transportation

The device may only be transported as described in section "Transportation of the device", page 104.

Operating position

The device may only be operated in the position shown in section "Preparation for commissioning", page 52.

Exchanging parts

Reliable functioning of the device and consequently safety for patients are only achieved if original spare parts from Kaltenbach & Voigt GmbH are used. Spare parts are manufactured by Kaltenbach & Voigt GmbH under especially high conditions with regard to material and manufacturing process and must be approved for use in a medial device in accordance with the relevant industrial standards.



Combination with accessories

The device may only be combined with approved fibre optic conductors and devices. For the selection of accessories, see chapter "Scope of delivery and accessories", page 30.

Sterilisation of accessories

The legal directives must be observed for the sterilisation of accessories.

Te instructions and warnings contained in the instructions for use of accessories must be observed.




Warning and instruction signs

Numerous warnings and instructions signs are attached to the device.

Figure 3 shows the warning and instruction plates and their position on the side and bottom of the device, Figure 4 shows the plates on the front of the device.

Figure 3. Warning and instruction plates (side and bottom of device)






No. in Figure 3	Sign		
1	SICHTBARE UND UNSICHTBARE LASERSTRAHLUNG BESTRAHLUNG VON AUGE ODER HAUT DURCH DIREKTE ODER STREUSTRAHLUNG VERMEIDEN LASER KLASSE 4 (IEC60825-1-2007-03)	Warning laser beam	
2	LASERKLASSE 4 LASERSTRAHLUNG DIODE GainAs 980±15nim SPEZIFIKATIONEN LASERLEISTUNG max 20W cw (nach IEC60825- 1:2007-03) PILOTLASER 640±10nm LASERLEISTUNG max. 1mW	Specification laser beam	
3	1	Foot switch connection	
4		Connection for door contact	
(5)	100 -240 VAC	Connection for power cable	
5	24 VDC	Connection for battery charger	
Kaltenbach & Voigt GmbH Bismarckring 39 D-88400 Biberach/Riß Type: MASTER lase		Type of device, For current manufacturing date, serial number and performance / supply data see nameplate on device	
6	2012-05 Ser.No.: D12M-0100 24VDC 60W Class: 2 Typ: BF 🕅 IPX1 Kaltenbach & Voigt GmbH Bismarckring 39 D-88400 Biberach/Riß	See "Product service life and disposal", page 111.	
		Device of protection class II according to EN 60601-1.	
		Type BF according to EN 60601-1.	



Figure 4. Warning and instruction plates, front of device



No. in Figure 4	Plate / symbol	
1		Follow the instructions for use!
2		Caution laser beam!
3		Fibre optic conductor connection



Description

Technical specifications

Model designation

RAVO EXPERT TASE	for use in the dental field
Electrical connection data	
Voltage	100 - 240 VAC
Frequency	50/60 Hz

KoVa EVDEDT loca I again devices with internal newer comply

Power consumption ≤ 120 VA

Power consumption:(max.) 2 A Climate (operating state)

Ambient temperature +15° C to +30° C

Relative humidity 30 % to 85 %

Air pressure 800 hPa to 1060 hPa

ON-time

P < 4 Watt: 66%, max. 15 min cw P = 4 - 5 Watt: 30%, max. 3 min cw

Dimensions and weight

	Cupboard	with display
Height	70 mm	85 mm
Width	142 mm	142 mm
Depth	145 mm	200 mm
Weight	approx. 1.4 kg	

Noise emission

Maximum $L \le 60 \text{ dB}$ in standby L < 50 dB

Fibres

Min. core diameter 240 μm na. ≥ 0.22

Therapy laser

Wavelengths 980 ± 15 nm

Laser power 0.2 - 8 Watt (see also Setting ranges, page 44)

Pilot laser

Wavelengths $640 \pm 10 \text{ nm}$ Laser power 0 to 1 mW



Model designation

KaVo MASTER lase	Laser devices with internal battery for use in the dental field
Electrical connection data	Charger adaptor

Power consumption: (max.) 3 A

Voltage 100 - 240 VAC, 24 VDC

Frequency 50/60 Hz Power consumption ≤ 300 VA

Climate (operating state)

Ambient temperature +15° C to +30° C Relative humidity 30 % to 85 % Air pressure 800 hPa to 1060 hPa

ON-time

P < 4 Watt: 66%, max. 15 min cw P = 4 - 6 Watt: 30%, max. 3 min cw P = 6 - 10 Watt: 20%, max. 2 min cw

Battery

Laser operation 5 W > 30 min.; 8 W > 15 min.

Standby (display switched on) > 3 hours Charging duration ≤ 30 min. (90%)

Dimensions and weight

	Cupboard	with display
Height	70 mm	85 mm
Width	142 mm	142 mm
Depth	145 mm	200 mm
Weight	approx. 1.4 kg	

Noise emission

Maximum	L ≤ 60 dB
in standby	L < 50 dB

Fibres

Min. core diameter 240 µm na. ≥ 0.22

Therapy laser

Wavelengths 980 ± 15 nm Laser power 0.2 - 10 Watt (see also Setting ranges, page 44)

Pilot laser

Wavelengths 640 ± 10 nm Laser power 0 to 1 mW

Accuracy of the defined values is compliant with EN 60601-2-22. Further technical details can be found in the current service manual for MASTER/EXPERT lase.



Scope of delivery and accessories

Different device variants can be selected. The accessories specified for the device are dependent on the model (see "Combination with accessories, page 24).

Figure 5. Arrangement of the transport case with instructions for use and accessories in a transport carton




No. in Figure 5	
1	Instructions for use
2	Bare fibre tips
3	Transport case containing device and accessories



⚠ CAUTION

The device case may only be opened in the position shown in Figure 6. Observe the upper / lower alignment of the corporate lettering!

Figure 6. Opening the device case





Scana of dolivory	Device variants		
Scope of delivery	=\\D=D=\	MASTER lase	
	EXPERT lase	Standard	Wireless
with internal power supply	1 x	-	-
with internal battery	-	1 x	1 x
Charger adaptor*	-	1 x	1 x
Wired foot switch	1 x	1 x	-
Wireless foot switch	-	-	1 x
Tool** for the wireless foot switch	-	-	1 x
Bare fibre tip 300 μm	3 x	3 x	3 x
Bare fibre tip 240 μm	2 x	2 x	2 x
Sapphire tip	-	2 x	2 x
Cleaning cap for sapphire tip	-	1 x	1 x
Tool for sapphire tip	-	2 x	2 x
Handpiece for bare fibre tip	1 x	1 x	1 x
Handpiece for sapphire tip	-	1 x	1 x
Transmission fibre optic conductor	1 x	1 x	1 x
Safety goggles (patient)	1 x	1 x	1 x
Safety goggles (operator)	2 x	2 x	2 x
Device case	1 x	1 x	1 x
Laser warning label	1 x	1 x	1 x

Article numbers for device accessories

Name	piece	Art. no.
Charger adaptor	1	1.009.5371
Safety goggles (patient)	1	1.009.5369
Safety goggles (operator)	2	1.009.5370

Article numbers for application accessories

Name	piece	Art. no.
Sapphire tip (2x) including cleaning cap (1x)	1	1.009.5368
Tool for sapphire tip	1	1.009.5365
Handpiece for sapphire tip	1	1.009.5363
Bare fibre tip 300 μm	20	1.009.5366
Bare fibre tip 240 μm	20	1.009.5367
Handpiece for bare fibre tip	1	1.009.5360
Transmission fibre optic conductor MASTER lase	1	1.009.5357
Transmission fibre optics conductor EXPERT lase	1	1.009.5359

^{*} The charging adaptor charges the battery (version with battery) and charges the battery of the wireless foot switch.

^{**} The tool is required for the assembly of the wireless foot switch for the initial commissioning (see Assemble wireless foot switch, page 55).



Figure 7. Arrangement of the EXPERT lase with accessories in the transport case



Placeholder for instructions on delivery







Figure 8. Arrangement of the MASTER lase with accessories in the transport case



Placeholder for instructions on delivery







Figure 9. Arrangement of the MASTER lase (option: wireless foot switch) with accessories in the transport case



Placeholder for instructions on delivery







Setup

The device is enclosed in portable housing.



- 1 Laser Stop: Emergency switch off of device
- 2 Switch on the device, access to menu settings

Switching on the device,

- Switching over the device:
 - Ready → Standby → Ready (for version with wired foot switch);
 Ready → Standby (for version with wireless foot switch);
- 4 Contact touchscreen (hereinafter referred to as touchscreen)
- (5) Connection for transmission fibre optics
- 6 Transmission fibre optic conductor winding (fibre winding)
- $\left(egin{array}{c} 7 \end{array}
 ight)$ Handpiece or fibre holder
- 8 Wired foot switch
- 9 Wireless foot switch wireless foot switch (only for device versions with wireless foot switch)
- Charging adaptor for charging the battery (version with battery) and for charging the battery of the wireless foot switch

Control panel




Terminals

Electrical terminals

Electrical terminals are located on the rear side of the device.

Figure 10. Connections (rear side)



Transmission optic fibres connection


The fibre optic conductor terminal is located on the side of the device and is for the connection of the transmission of the transmission fibre optic conductor.

Figure 11. Device with transmission fibre optic conductor





Operating elements of the wireless foot switch



No.	Name	Function	
1	Left button	Switches the device from "Standby" mode to "Ready" mode and back.	
2	LEDs	LEDs shining indicates that the device is in "Ready" mode	
3	Middle button	Triggers the therapy laser beam.	
4	Right button	Switches the performance of the device one level higher	



The wireless foot switch id equipped with a position sensor. If the wireless foot switch is considerably inclined for longer than 3 sec., the device emits an error message. If this position is not rectified within 5 sec., the foot switch logs off automatically. You must log the wireless foot switch into the device again to continue the therapy (see Assemble wireless foot switch and log into the device, page 55).



Power supply of the wireless foot switch

The supply to the wireless foot switch is secured by a battery. The scope of delivery includes a charger adaptor (Scope of delivery and accessories, see page 30).



If the battery charging is inadequate, a symbol will appear at the bottom of the screen on the left. From this time, the foot switch will remain in operation for approx. 5 hours. Charge the battery after completing the running therapy.



You can also read the status of the battery on the screen via the menu: menu \rightarrow wireless foot switch \rightarrow status.



To charge the battery, connect the wireless foot switch to the power supply by means of the charger adaptor included in the delivery.



When connecting the charger adaptor with the foot switch, check the position of the groove and spring of the plug connection (red mark on the plug must be pointing upwards).



The charge symbol is blended out when the battery has been adequately charged.

The course of therapy is not disturbed by this procedure.



Power supply of the battery

The scope of delivery of the battery-operated device includes a charger adaptor (Scope of delivery and accessories see page 30).



If the battery charging is below 50%, a symbol will appear at the bottom of the screen on the right.

For information concerning the remaining duration of operation, see chapter "Technical specifications", page 28.



To charge the battery, connect the device (see Terminals, page 37) to the power supply by means of the charger adaptor.

The charging process is indicated by the arrow in the symbol.

For information concerning the remaining duration of operation, see chapter "Technical specifications", page 28.

The course of therapy is not disturbed by this procedure.



To extend the duration of operation, chose a brief automatic switch-off time (see "Setting the switch-off time", page 89) or switch the device off as soon as it is no longer needed.



Control panel

The control panel serves the input and display of the device functions (operating modes, setting parameters etc.) and consists of the display and the touch keys.

The display is equipped with a touch-sensitive surface (touchscreen). the touchscreen is divided into three areas which are emphasised by different background colours.

The input is made by pressing the keys and touching a displayed element or value.

⚠ WARNING

Undesired operating functions are possible! If several pressure points are activated, this could lead to an undesirable operating function.

On input, it must be ensured that only one pressure point on the touchscreen is active.

⚠ CAUTION

Make entries by gently touching the touchscreen, excessive pressure can lead to destruction of the display.

Figure 12. Control panel





No. in Figure 12	Name	Function
1	Lower status bar	The device status is displayed in the status bar with messages. A change of status leads to an immediate update of the status display. The status bar also serves the display of messages as well as the release of the foot switch for laser activation. Activation of the status bar releases Ithe fioit switch to trigger the laser 2 seconds after activation.
2	Slider scale	Setting the indiviudal laser parameters. The activated value can be changed by direct activation of one of the values displayed in the slider scale. For confirmation of the activation, the colour of the selected element or value changes from white to green. The range of values available for selection can be shifted to the left or right or displayed by "dragging" over the slider scale.
3	Laser parameter field	Laser parameters are displayed and can be changed by direct activation of one of the values displayed in the slider scale. The slider scale is blended in on selection of a corresponding laser parameter.
4	Upper status bar	The therapy status is displayed in the status bar. The therapy status can be changed via the drop-down menu. The respective drop-down menu is opened on selection of a corresponding status element.
5	Indication element	Setting, call and programming of the der application oriented laser parameters.
6	Pulse form element	Selection of the pulse form of the therapeutic laser beam.
7	Protocol element	Selection of the displayed protocol data and timers as well as resetting all protocol data and timers to 0.
8	Pilot element	Setting the pilot brightness level.
9	Switch-on button	Switching on the device, Switch over Ready < - > Standby (see page 36).
10	Laser Stop button	Emergency device switch off
11	Menu button	Switch on the device, access to menu settings.



Therapy mode Setting parameters

Table 1. Adjustable pulse forms

Pulse forms		Description		
Continual: (cw)		The set laser performance is emitted as soon as the foot switch is activated.		
Continual with setting of the application time:		The set laser performance is emitted for a pre- selected time as soon as the foot switch is activated.		
Interval:		The set laser performance is emitted in a pre- selected pulse pause / pulse duration ratio as soon as the foot switch is activated.		
Interval with setting of the application time:	Me	The set laser performance is emitted in a pre- selected pulse pause / pulse duration ratio for a pre-selected time as soon as the foot switch is activated.		
Peak:		The maximum performance of the laser device (Pp) is emitted in a pulse pause / pulse duration ratio as soon as the foot switch is activated. Pulse duration (tp) and average performance (Pa) are set. The necessary pulse pause (ti) is calculated automatically by the system from the selected pulse duration (tp), the average performance (Pa) and the maximum performance (Pp) in µs accuracy. e.g.: tp=1ms, Pa=5W, Pp=10W → ti = 1ms		

Table 2. Adjustable parameters

Parameters	Description
Power	Setting the power for the therapy laser beam.
Pulse duration	Setting the active phase of a pulse of the therapy laser beam.
Pulse pause	Setting the inactive phase between two pulses of the therapy laser beam.
Duration	Setting the application time for the therapy laser beam.



Table 3. Setting ranges

Tallot of Cottaining Talling Co						
		Laser power				
Device variants	Parameters			M		
EXPERT	Power	0.2 - 5W	0.2 - 5W	0.2 - 5W	0.2 - 5W	0.2 - 4W*
lase	Peak performance	-	-	-	-	8W**
MASTER	Power	0.2 - 8W	0.2 - 8W	0.2 - 10W	0.2 - 10W	0.2 - 5W*
lase	Peak performance	-	-	-	-	10W**
		Duration				
EXPERT/	Pulse duration	cw	cw	0.5ms - 3s	0.5ms - 3s	0.5 – 50ms
MASTER lase	Pulse pause	cw	cw	0.5ms - 3s	0.5ms - 3s	-
	Timer	-	0.1 – 90s	-	0.1 – 90s	-

^{*} adjustable average power

^{**} output peak power



Selection of indication

The treatment parameters are specified by selection of the desired indication (recommendations).

Treatment instructions compiled by experienced specialists and reference physicians are given with the indications. Special cautionary measures are prefixed with CAVE and displayed in orange.

⚠ WARNING

All information is kept as conservative as possible and has been provided with the utmost care and according to our best knowledge and judgement. The achievement of the required effect on tissue is to be obtained by adaptation of the laser parameters if necessary and must, in any case, be verified by suitable measures. The treatment instructions of the current specialised literature should be consulted for this purpose. Kaltenbach & Voigt GmbH shall assume no liability for personal injury and materiaL damage resulting from application of indication settings and indication instructions.

⚠ WARNING

Irreversible damage to health could be caused in soft and hard tissue. MASTER laser or EXPERT lase must be released for use with the bleaching material by the manufacturer of the bleaching material. Comply with the instructions for use of the manufacturer of the bleaching material.

⚠ WARNING

Irreversible damage to health could be caused in soft and hard tissue.

Always keep the laser fibre in motion and do not concentrate one one point. If possible, always laser parallel to the tooth and periodontium.

The indications can be saved in a list of Favorites and the settings changed to individual requirements (see chapter "Compilation and use of list of Favorites", page 81).

The performance is restricted in some specialist areas for safety areas.



Table 4. Selection of indications

Table 4. Selection of indications						
	Maximum		Indication			
Specialised areas adjustable performance		Display on indication selection	Display on the indication element button	Full text		
Surgery						
		Abscess	Abscess			
		Adenoma	Adenoma			
		Biopsy	Biopsy			
		Drainage	Drainage			
		Epulis	Epulis			
		Fibroma	Fibroma			
		Fistula duct	Fistula duct			
		Labial frenectomy	Labial frenectomy			
		Lingual frenect.	Lingual frenect.			
	8 W	Gingivectomy	Gingivectomy			
		Gingivoplasty	Gingivoplasty			
		Hemangioma	Hemangioma			
		Hemostasis	Hemostasis			
		Hyperplasias	Hyperplasias			
		Implant exposure	Implant exposure			
		Incision / excision	Incision / excision			
		Leukoplakias	Leukoplakias			
		Lichen planus	Lichen planus			
		Mucocele /ranula	Mucocele /ranula			
		Operculectomy	Operculectomy			
		Papilloma	Papilloma			
		Vestibuloplasty	Vestibuloplasty			



	Maximum adjustable performance	Indication			
Specialised areas		Display on indication selection	Display on the indication element button	Full text	
Endo / Paro					
		Endodontics	Endodontics	Endodontics adjuvant	
		Endo periapical	Endo periapical	Endo adjuvant in periapical defect	
		Exstirpation	Exstirpation		
	4 W	Periodontology	Periodontology	Periodontology adjuvant	
		Periimplantitis	Periimplantitis	Peri-implantitis decontamination of the implant	
		Pulpotomy	Pulpotomy		
		Pulp capping	Pulp capping		
Others					
		Aphthea, herpes	Aphthea, herpes		
		Bleaching	Bleaching	Activation of bleaching material by means of laser light	
	8 W	Desensitisation	Desensitisation		
		Mycoses	Mycoses		
		Sulcus dilation	Sulcus dilation	Low Lovel	
		LLLT	LLLT	Low Level Laser Therapy	
		Blacken the fiber	Blacken the fiber		
		Customized 1 to 5	Call indications		



Pilot

The pilot laser simulates an invisible therapy laser beam by means of visible light. The brightness of the pilot laser can be set from 0 to 5 in single increments.

The pilot laser is a laser of the laser class 2 according to IEC EN 60825:2007-03.



If the pilot is switched off and invisible, the device may only be used for contact applications.

Protocol

The device logs and displays automatically:

- the sum of the applied energy and application times
- the duration of the inactive phase of the laser beam
- the duration of the inactive phase and the duration of the active phase of the laser beam
- Peak performance and average performance

You will find information concerning the selection of the Changing the display mode of the protocol button on page 86.

The protocol data can be reset to 0 via the protocol element.



The protocol data are principally reset to 0 when the device is switched on.



Menu

.The menu enables the individual definition and saving of basic settings for the user.

Table 5. Adjustable parameters in the menu

Fun	ction / setting	Description
- Se	ettings Tab	
	/olume	Volume setting
	Laser audible	Volume of the audible warning
	warning	the volume of the audible warning can be set at 5 levels.
	Input feedback	volume of the feedback tone on button activation the volume of the audible warning can be set at 5 levels.
7	Tool Option	
	Switch-off time (min.)	Set the time for the automatic switch off. Devices with mains power supply: 3, 5, 10, 15, 30 minutes or Off. Devices with battery supply: 3, 5, 10, 15 minutes. The devices switches off automatically if it is not used
		over a time period.
	Transitional period Ready → Standby (min.)	The time required for setting the automatic status change of the device:1, 2, 3 minutes.
		The devices switches automatically to STANDBY mode if it is not used over a set time period when it is the READY mode.
		Switching the pilot laser on and off
	Pilot during treatment	If you have set the the selection <i>Off</i> for the <i>Pilot</i> option <i>during the treatment</i> , the pilot laser is switched off automatically when you activate the foot switch.
C	Operation	Select operating mode
	Therapy	Treatment is possible in this setting
	Demo	Demonstration of functions and operation only.
	Service	Only for authorised service personnel
	Switch off	Turn device off



	Language	Select language for display texts			
	Deutsch	Texts in the display appear in German			
	English	Texts in the display appear in English			
	Français	Texts in the display appear in French			
	Español	Texts in the display appear in Spanish			
	Italiano	Texts in the display appear in Italian			
	Portuguese	Texts in the display appear in Portuguese			
	Svenska	Texts in the display appear in Swedish			
	Norsk	Texts in the display appear in Norwegian			
	Polski* Texts in the display appear in Polish				
	Time	Set the real-time clock			
	Clask	Time (hour/minute)			
	Clock	The time is specified in 24/hour format.			
	Data	Date (day/month/year)			
	Date	The date is specified in dd.mm.yy format.			
	Pin	Setting the frequency of the automatic PIN - request (Off, daily, ON)			
		Change the access code			
Lo	gbook				
	Up to 255 treatments can be recorded in the logbook. A new entry is made when the data changes and laser performance is applied.				
	Operation	Read operating logbook			
	Therapy Reading therapy logbook				
W	ireless foot switch				
	Login	Instructions for logging the foot switch into the device			
	Status	Information about the foot switch status			

^{*} from SW - Version 1.3.1



Operation

This chapter assumes that you are familiar with the operating elements of the device.



Read and observe the chapters "Safety" "Description".

Observe safety

Carry out a visual inspection prior to every startup:

⚠ CAUTION

You may not start the device if the safety checks give reason for concern.

In this case, notify your responsible after sales service.

- 1. Check the control panel and the housing of the device for mechanical damage.
- **2.** Check the existing cable for mechanical damage.
- **3.** Check the existing connections for mechanical damage.



Preparation for commissioning

A wide range of devices and accessories are available for the MASTER/EXPERT lase.



Observe the instructions for use of the devices.

⚠ WARNING

If the laser beam hits inner parts of the device, this could damage the eyes of the user and damage the device. The fibre optic conductor must protrude distally on laser activation from the device canal.

⚠ CAUTION

The MASTER/EXPERT lase is an device of protection class IPX1 in accordance with IEC60529, i.e. it is equipped with normal protective measures against the ingress of dust and fluids and is not waterproof.

Malfunctioning could occur in case of high air humidity or ingress of water and dust into the housing. Corrosion of the internal mechanical and electronic systems can lead to irreparable damage.

Condensation can from on the optic components as a result of sudden changes in temperature (e.g. relocation of the device from a warm area with high air humidity to an air conditioned cooler room).

It is essential to ensure, prior to startup, that the MASTER/EXPERT lase is set up several hours in the intended operating area in order to prevent destruction of optic components due to condensation.

⚠ CAUTION

On application of the device with battery:

Charge the battery of the device completely (approx. 1 hour). Observe the information in chapter "Power supply of the battery", page 40).

On application of the wireless foot switch:

Charge the battery of the foot switch completely (approx. 5 hours) prior to the initial commissioning. Observe the information in chapter "Power supply of the wireless foot switch", page 39).



Only use the charger adaptor provided to charge the battery.






Mounting the charger adaptor









Dismounting the charger adaptor







⚠ CAUTION



Observe the position of the plug when you connect the charger adaptor:

The marking (red dot) on the sleeve of the plug must be pointing downwards when you connect the device.

The marking (red dot) on the sleeve of the plug must be pointing upwards when you connect the device to the wireless foot switch.




1. Bring the device into the operating position as shown in Figure 13.

Figure 13. Operating position



- 2. Clearly define and mark the laser area, see section "Laser protection officer".
- Only for models with mains power supply: Connect the power cable to the power supply.
- 4. Connect the foot switch (wired) or log the wireless foot switch into the device.

⚠ CAUTION





Observe the position of the plug.

The marking (groove) on the sleeve of the

The marking (groove) on the sleeve of the plug must be pointing upwards when you connect the footswitch to the device.




Assemble wireless foot switch and log into the device, page 55.




Assemble wireless foot switch and log into the device



The steps 1 to 3 must be carried out only for the initial commissioning of the wireless foot switch. When the wireless foot switch is assembled, begin with step 4.



1. Remove the wireless foot switch parts from the packaging.



Guide the ends of the bracket into the intended recesses at the bottom of the wireless foot switch.



Screw both ends of the bracket tightly to the foot switch, fastening the sealing cap of the bush.



The screws and required tolls are included in the scope of delivery (Scope of delivery and accessories, page 30).

Store the tools in the devices case after the assembly.



4. Place the wireless foot switch near to the device.





The wireless foot switch may not be connected to the device and the device must be switched on. When switching on the device, observe chapter "Commissioning", page 57.



5. Activate the lower status bar.



6. Follow the instructions shown on the screen.



7. Confirm selection with





The device changes to therapy mode. The symbol for the wireless foot switch is blended in below, on the left.




if you possess several KaVo devices with wireless foot switches, you have to log the foot switch into the device each time you switch on.



Commissioning




Use of the operating units or setting possibilities in any other manner as described here can lead to dangerous radiation.



 Carry out the preparation for startup.

Open the protective bracket of the foot switch (by gently touching the bracket).



switch the device on



KaVo MASTER lase Dental Diode Laser 980 nm, 10 W SW Version X.X.X (XXX)

For approx. 3 seconds, the startup image is displayed at the left of the screen (example: EXPERT lase), the devices performs a self test.

The number in brackets is only relevant for internal administration purposes.



if a safety check is required, this switchon image appears for approx. 3 seconds (example).

Warning notes and instructions for possible maintenance are displayed.





if an error occurs during the startup, the device will show this image with an ERROR number and an explanatory text. In this case carry out a new start. If the error persists, please notify the KaVo Service and tell the service employee the error message, the software version and the serial number of the device.



If the self test has been completed successfully, the screen for the entry of the access code (PIN) on the device is blended in.

The cursor is displayed by a green dot.



With number input the green dot is replaced by *.

Enter the PIN via the number bar (the PIN is a default setting 1234).

The number bar is blended out if the identification is successful.



After approx. 1 sec., this screen is blended in on the device (the display of the screen is an example).



Every laser device in current disuse must be safeguarded against unauthorised use.

Changing the factory PIN (see Changing the PIN, page 96).

if you have forgotten your personal PIN, you can start up the device with a Super PIN. You can obtain the Super PIN from Kaltenbach & Voigt GmbH Service Support. After the start with the Super PIN, change the PIN on your personal PIN again.



Commissioning of the transmission fibre optic conductor


⚠ WARNING

Follow the instructions for commissioning and sterilisation in the Instructions for use of the relevant accessory.

⚠ WARNING

Danger of contamination!

Carry out the wipe disinfection of the transmission fibre optic conductor both just before application as well as directly after application including the initial commissioning.



Connect the transmission optic fibres to the laser device.



The plug of the fibre optic conductor must be screwed on completely, approx. 5 clockwise turns are required.

Verify the firm seating of the fibre optic conductor plug.

⚠ CAUTION

To prevent damage to the fibre optic conductor:

hold the transmission fibre optic conductor firmly in the guide grove while winding (a);

while winding, grip the transmission fibre optic conductor at a distance of 1 cm to 2 cm from the edge of the fibre winding (b, c), so that the fibre optic conductor can be wound gently under the white silicone disc;







Avoid contact with the distal and proximal ends of the transmission fibre optic conductor.

The part of the transmission fibre optic conductor, which is not required for the treatment, can be completely wound under the fibre winding (bending radius is 25 mm).





 Wind or unwind the transmission fibre optic conductor on the fibre winding according to requirement, ensuring that the length of the conductor is adequate for the treatment.



Put the housing cover on, guiding the transmission fibre optic conductor through one of the recesses provided for this purpose.







- 3. Check the eye protection.
- 4. Switch on the laser device.

The laser device detects the connection of the fibre optic conductor.

The message "Connect the fibre optic conductor to the device" extinguishes. The laser status bar is blended in.



⚠ WARNING


Cross contamination hazard!

Do not put on the protective cap during applications.

The protective cap serves as protection only for delivery or in case of prolonged storage in the case.

5. If necessary, remove the protective cap of the fibre optic conductor.






6. Direct the distal end of the transmission fibre optic conductor (distance: approx. 20 mm) vertically to a bright, flat surface and observe the light spot of the pilot laser beam (radiation characteristics).

In case of a technically faultless transmission fibre optic conductor, an even circular light spot is generated with intensity rapidly decreasing to the outside.

⚠ WARNING

If the radiation characteristics do not comply with the above mentioned description, or presentation, the complete functioning of the fibre optic conductor tip cannot be achieved. Hence the desired success of treatment cannot be guaranteed or the handpiece is overheated. If necessary, replace the fibre optic conductor and check again.

Figure 14. Device with connected transmission optic fibre conductor and div. accessories





Mounting the sapphire tip handpiece and commissioning



The handpieces and sapphire tips are delivered in non-sterile condition by the manufacturer.



To enable sterile application, handpieces and sapphire tips must prepared prior to every application including the initial commissioning together with the tools (see Instructions for Use Accessories for Dental Application).

⚠ WARNING

It is essential to prevent contamination of the sterilised handpieces and sapphire tips that come into contact with patients.

Observe the instructions for the aseptic removal and assembly of the handpieces and sapphire tips.



Contamination hazard!

If sterile packaging is damaged, previously sterilised handpieces and sapphire tips must be sterilised again before use.



Contamination hazard!



If possible, keep the handpiece in your hand and do not lay it down on an unsterile surface. Avoid contact with the sapphire tip.


Use only a sterile holder throughout the entire process!






Transmission optic fibre conductors must be connected to the device by means of fixing sleeves (see page 59)

1. Open the pack at the handle side of the handpiece and remove the fixing sleeve (1a) and the chuck (1b).



Guide the transmission optic fibre conductor through the opening in the fixing sleeve.



3. Guide the optic fibre conductor through the opening in the chuck and screw the chuck onto the fixing sleeve (2) of the optic fibre conductor.



The chuck must be hand-tightened onto the fixing sleeve of the optic fibre conductor. Verify the firm seating of the chuck before you continue.



Mount the fixing sleeve, applying slight pressure to the chuck.





5. Remove the handpiece from the packaging.



6. Open the sapphire tip pack at the tool side.



Remove the tool with sapphire tip from the pack, avoiding contact with the sapphire tip.



8. Screw the sapphire tip into the handpiece in a clockwise direction.

9. Pull off the tool from the sapphire tip.









10. Guide the transmission fibre optic conductor into the handpiece to backstop and screw the fixing sleeve to the proximal end of the handpiece tightly.



Guide the transmission fibre optic conductor slowly through the curve of the handpiece, exercising only slight pressure. Excessive pressure compresses the transmission optic fibre conductor and intensifies the resistance.

- 11. Check that there is no gap at the screw connection between the clamping and the handpiece and verify the firm sea of the transmission optic fibre conductor.
- 12. Switch the laser device on and verify that the pilot light an the distal end of the sapphire tip is clearly detectable. If necessary, check the screw connections of the clamping on the fixing sleeve and on the handpiece.
- **13.**Set the desired parameters on the laser device and commence the treatment.

In case of a technically faultless sapphire tip, an even light spot is generated with intensity rapidly decreasing to the outside.

⚠ WARNING

If the radiation characteristics do not comply with the above mentioned description, or presentation, the handpiece could heat up inadmissibly and the complete functioning of the sapphire tip cannot be achieved. Hence the desired success of treatment cannot be guaranteed. If possible, keep the handpiece in your hand and do not lay it down on an unsterile surface.



Mounting the bare fibre handpiece and commissioning

⚠ WARNING

The bare fibre handpieces are delivered in non-sterile condition by the manufacturer.

To enable sterile application, bare fibre handpieces must prepared prior to every application including the initial commissioning (see Instructions for Use *Accessories for Dental Application*).

⚠ WARNING

It is essential to prevent contamination of the sterilised handpieces and bare fibre tips that come into contact with patients.

Observe the instructions for the aseptic removal and assembly of the handpieces and bare fibre tips.

⚠ WARNING

If sterile packaging is damaged, previously sterilised handpieces must be sterilised again before use.

⚠ WARNING

Contamination hazard!

If possible, keep the handpiece in your hand and do not lay it down on an unsterile surface. Avoid contact with the bare fibre tip.

Use only a sterile holder throughout the entire process!


STERILE EO



The bare fibre tips are delivered in sterile condition by the manufacturer (EO = ethylene oxide).



The bare fibre tips from damaged sterile packaging may not be used.




the transmission optic fibre conductor must be connected to the device (see page 59)

 Open the pack at the handle side of the handpiece and remove the fixing sleeve (1a) and the chuck (1b).













2. Mount the fixing sleeve, applying slight pressure to the chuck.

- 3. Guide the transmission optic fibre conductor through the opening in the fixing sleeve and the chuck.
- 4. Open the pack at the adaptor side of the bare fibre tip.
- 5. Remove the bare fibre tip from the pack, avoiding contact with the fibre of the fibre optic conductor.



Hold the bare fibre tip without touching the snap-in hook.

- 6. Attach the adaptor of the bare fibre tip to the distal end of the transmission fibre optic conductor until the adaptor perceptibly snaps into place.
- 7. Remove the handpiece from the pack, avoiding contact with the distal end of the handpiece.





8. Guide the transmission fibre optic conductor with the mounted bare fibre tip into the handpiece.



Screw the fixing sleeve at the proximal end of the handpiece so that the fibre optic conductor will still have freedom of movement.




10. Adjust the required length of the bare fibre tip at the distal end. The distal end of the fibre optic conductor must protrude from the handpiece by at least 5 mm.




11. Screw the clamping at the proximal end of the handpiece tightly.





12.Verify the firm seating of the transmission fibre optic conductor.



13. Switch on the laser device. Direct the distal end of the bare fibre tip (distance: approx. 20 mm) vertically to a bright, flat surface and observe the light spot of the pilot laser beam (radiation characteristics). If necessary, check the screw connections of the clamping on the handpiece.

In case of a technically faultless bare fibre tip, an even circular light spot is generated with intensity rapidly decreasing to the outside.

⚠ WARNING

If the radiation characteristics do not correspond to the abovementioned description or depiction, the handpiece might heat up inadmissibly, which may lead to burn injury. You may need to replace the BareFiber-Tip and check the new tip again.

⚠ WARNING

Undesirable side effects could arise in case of contact cutting (e.g. overheating of deeper layers of tissue).

To prevent complications, the sapphire tip or bare fibre tip must be burnt in (see Blacken the sapphire tip / bare fibre tip, page 70).

14.If you have selected an indication that is not performed in the contact, set the desired parameters on the laser device and commence the treatment.



Blacken the sapphire tip / bare fibre tip



1. Activate the indication element.



2. Activate \bigcirc the display others.

The display *Others* changes colour to green and moves to the left.



 "Draw" up via the indication list and select Blacken the fiber from the list.



4. Confirm with OK to proceed to therapy mode.



The parameters for *Blacken the fiber* are set automatically.







5. Touch the imprinted dark square on the sapphire tip & bare fibre tip pack with the distal end of the sapphire tip / bare fibre tip and activate the laser beam.





6. Check that the sapphire tip / bare fibre tip has been burnt in (a dark coating must be detectable at the distal end of the sapphire tip / bare fibre tip), if necessary, repeat the procedure from step 5.

- **7.** Set the desired indication on the laser device.
- 8. Commence the treatment.



Dismounting the sapphire tip handpiece



 Unscrew the fixing sleeve from the handpiece completely and pull the transmission fibre optic conductor with the chuck and fixing sleeve slowly from the handpiece.



2. Unscrew the chuck from the fixing sleeve and pull the chuck with the fixing sleeve from the transmission fibre optic conductor.



3. Unscrew the sapphire tip from the handpiece using the tool.



4. Pull off the sapphire tip.

 Then carry put the preparation process (see Instructions for use Accessories for dental application).



Dismantling the bare fibre handpiece




1. Unscrew the fixing sleeve from the bare fibre handpiece and pull the transmission fibre optic conductor with the bare fibre tip from the handpiece.





Hold the bare fibre tip without touching the snap-in hook.

2. Pull the adaptor of the bare fibre tip from the distal end of the transmission fibre optic conductor.



- 3. Pull the chuck with the fixing sleeve from the transmission fibre optic conductor.
- 4. Then carry put the preparation process (see Instructions for use Accessories for dental application).



At the end of the service life cycle, dispose of the optic fibres and handpieces in accordance with the local regulations for contaminated products. In this way you will minimise hazards for the environment and personnel that could be caused by contamination and residues of the fibres and handpieces employed.



Therapy settings

To be able to monitor the set parameters during therapy, secure a clear view of the display of the device.

Select the indication



Unsuitable fibre optic conductors could be damaged.

Only use suitable fibre optic conductors when you change the indication or the specialist area. Follow the instructions in the Instructions for use of the fibre optic conductor.



1. Activate the indication element e.g. Abscess.



Select the specialist area (e.g. Endo/Paro).

The display for the chosen specialised area changes colour to green and moves to the left. The corresponding indications are blended in.



Select the indication (e.g. Pulpotomy).

White arrows are blended in if the list of indications is longer than can be shown on the screen. The indications can be shifted by "dragging" up or down over the indication list.



- 4. Read the instructions with care.
- **5.** Confirm selection with $^{\circlearrowleft}$ *OK* to proceed to therapy mode.

On return to therapy operation, the chosen indication is blended into the upper status bar and the parameters are set to the chosen indication.



Setting the performance



The desired indication must be chosen (Select the indication, page 74).



1. Activate the performance value.



The chosen value changes the colour from white to green and the slider scale is blended in.

2. Activate the desired value on the slider scale.



The chosen value is set.

if you want to set a value that is higher or lower than is indicated on the display, shift the slider scale in the corresponding direction.



After approx. 5s, the colour of the display changes from green to white and the slider scale is blended out.



Setting the performance with the wireless foot switch



The desired indication must be chosen (Select the indication, page 74).



1. Activate the right button on the wireless foot switch.

The chosen value changes the colour from white to green, the next highest value is set at once and the slider scale is blended in.

⚠ WARNING

Continually check the performance display visually!



2. Activate the right button on the wireless foot switch repeatedly until the desired value is activated on the slider scale.

The cursor (value marked in green) moves from left to right within the visible area of the slider scale.



When the highest value at the end of the visible area of the slider scale has been reached, the cursor springs to the first lowest vale at the beginning of the slider scale.



The chosen value is set at once.



After approx. 5s, the colour of the display changes from green to white and the slider scale is blended out.



Setting the pulse form



The desired indication must be chosen (Select the indication, page 74).



Activate the element Pulse form.



The activated element changes colour from white to green.

The pulse forms for selection (Drop-Down list) are blended in.

2. Select ** the desired pulse form in the Drop-Down list.



The symbol for the chosen (activated) pulse form changes colour to green.



The selected pulse form can be adopted actively by pressing the pulse form element. Otherwise the chosen pulse form is set automatically and the pulse forms for selection (Drop-Down list) are blended out after approx. 2 s. The element changes colour from green to white. The corresponding laser parameter values are displayed and can now be changed.



Setting the pulse duration













The desired indication must be chosen (Select the indication, page 74).

Confirm the value of pulse duration (provided with ... icon).

The chosen value changes the colour from white to green and the slider scale is blended in.

2. Activate the desired value on the slider scale.

The chosen value is set.

After approx. 5s, the colour of the display changes from green to white, the slider scale is blended out.



Setting the pulse pause











The desired indication must be chosen (Select the indication, page 74).

1. Confirm [→] the value of pulse duration (displayed with ¬ icon).

The chosen value changes the colour from white to green and the slider scale is blended in.

2. Activate the desired value on the slider scale.

The chosen value is set.

After approx. 5s, the colour of the display changes from green to white, the slider scale is blended out.



Setting the duration of application













The desired indication must be chosen (Select the indication, page 74).

1. Activate the value for the duration of application.

The chosen value changes the colour from white to green and the slider scale is blended in.

2. Activate the desired value on the slider scale.

The chosen value is set.

After approx. 5s, the colour of the display changes from green to white, the slider scale is blended out. Additional information can be obtained from: pages 43 and 44.



Compilation and use of list of Favorites Saving the default settings of the indication parameters



Each indication can be saved in the list of Favorites. The indication will be overwritten if it is already in the list of Favorites.



1. Activate ** the indication element.



2. Select the desired specialist area (** e.g. Endo/Paro*).

The display for the chosen specialised area changes colour to green and moves to the left. The corresponding indications are blended in.



3. Select the desired indication (e.g. *Periodontology*).



Activate

Save in Favorites.

The Favorites display changes colour to green and moves to the left. The desired indication is included in the list of Favorites.



 Confirm with SOK to proceed to therapy mode.



Saving individual indication parameter settings



Each indication can be saved with any settings (including pilot). The indication will be overwritten if it is already in the list of Favorites.



 Select the indication (see. "Select the indication", page 74),
 e.g. Paro adjuvant.

The recommended therapy parameters (default settings) are set.



2. Set the desired values.

As soon as you have changed the recommended therapy parameters, the colour of the indication element changes to yellow.

3. Activate ** the indication element.



- Activate
 - Save in Favorites.



The Favorites display changes colour to green and moves to the left. The changed indication is included in the list of Favorites.




5. Confirm with SOK to proceed to therapy mode.



Calling up indication parameters from the list of Favorites











- 1. Activate the indication element.
- 2. Activate Favorites.

The Favorites display changes colour to green and moves to the left.

3. Select the indication from the list of Favorites.

White arrows are blended in if the list of Favorites is longer than can be shown on the screen. The indications can be shifted by "dragging" up or down over the list of Favorites.

4. Confirm with SOK to proceed to therapy mode.

The indication parameters are set automatically.

the indication element is displayed in yellow if the indication has been saved with individual changed parameters.



Deleting indication parameters from the list of Favorites











1. Activate the indication element.

Select menu Favorites if necessary.

2. Select the indication to be deleted from the list of Favorites.

The name of the indication changes colour from white to green.

3. Activate Delete.

The indication is removed from the list of Favorites. The indication following in the list of Favorites is activated.

4. Confirm with $^{\circlearrowleft}$ OK to proceed to therapy mode.



Setting the brightness of the pilot laser



The fibre optic conductor must be connected to the device to be able to visually check the brightness of the pilot laser.





The desired indication must be chosen ("Select the indication", page 74).

1. Activate sthe element *Pilot*.



The activated element changes colour from white to green. The bar with the pilot brightness level is blended in.

The activated pilot brightness levels are displayed red.



2. Activate the pilot brightness levels on the bar that has been blended in.

The brightness setting can be checked visually on the connected fibre optic conductor.

The pilot laser can be switched off. To achieve this, activate the lower pilot brightness level twice.



Risk of injury!

If the pilot is switched off and invisible, the device may only be used for contact applications.



The selected value can be adopted actively by pressing the pilot element. Otherwise the chosen value is set automatically. After 2 seconds, the bar with the pilot brightness levels is blended and the colour of the element changes from green to white.



Changing the display mode of the protocol button



Pilot Fine 1.00

Pilot Fine 1.00

Energy, 194

Time 1.00

Energy, Time 9.00

Goff 9.00

Reset



1. Activate the protocol button.

2. Select ** the desired pulse lines in the Drop-Down menu.

Display Σ Energy, Time:

The sum of the applied energy and the sum of the application times are displayed.

Display 🖰 off (mm:ss):

the duration of the inactive phase of the laser beam is displayed.

Display \bigcirc off / on (mm:ss):

the duration of the inactive phase and the duration of the inactive phase and the duration of the inactive phase of the laser beam are indicated alternatively.

Display \hat{P}, \bar{P} :

The peak performance (\hat{P}) and average performance resulting from the parameter selection

 (\bar{P}) are displayed.



Resetting the protocol or timer



The data to be reset must be displayed in the upper status bar.



1. Activate The protocol element.



The activated element changes colour from white to green.



2. Activate **Reset* on the Dropdown menu blended in.

The saved data are deleted from the protocol. The protocol element changes colour from green to white. The drop-down menu is bended out.

When the device is switched off after completion of treatment, the saved data are deleted from the protocol automatically.



Menu settings



The setting in the menu is made by touching a displayed element or value. For confirmation, the colour of the selected element or value changes from white to green.

Setting the volume for laser warning signal tones and input feedback



Begin with step 2 if you are already in the settings menu.

Activate the "Menu" button.



2. Activate the laser warning signal / input feedback levels in the relevant bar.

The volume of the tone can be acoustically controlled.

Activate the lower volume level twice to switch off the input feedback.



3. Confirm the setting with



Setting the switch-off time



Volume
Ontions

Volume
Ontions

2 ration
Language
Time
PIN

Pilot during therapy
On

OK

RF foot-switch
RF foot-switch
RF foot-switch
RF foot-switch
RF foot-switch
RF foot-switch
Size A control of time (min.)
1 5 30
1 2 3

Off

OK





Begin with step 3 if you are already in the menu - settings \rightarrow options.

1. Activate the "Menu" button.

- **2.** Activate the menu dialogue *Options*.
- Select the desired switch-off time and activate the corresponding symbol (e.g. 515).

Here you can also deactivate the automatic switchoff of the device in device variants with an internal power supply (Off).

This possibility is not available in battery-operated devices.

4. Confirm the setting with $\mathcal{S} \cap \mathcal{K}$



Setting transitional period Ready → Standby



Begin with step 3 if you are already in the menu - settings \rightarrow options.

1. Activate the "Menu" button.



Activate the menu dialogue Options.



The time required for the automatic status change of the device: 1. 2. 3 minutes



Select the desired transitional period (e.g. 3) and activate the corresponding symbol.



4. Confirm the setting with \mathcal{S} *OK*.



Selecting option "Pilot during treatment"



Begin with step 3 if you are already in the menu - settings \rightarrow options.

1. Activate the "Menu" button.



2. Activate the menu dialogue
 Options.



3. Select the desired setting under the line "Pilot during treatment", to achieve this, activate the corresponding selection (e.g.
Off).



4. Confirm the input with

☞ OK.

the device changes to therapy mode (observe the information in chapter Pilo, page 48).



Selecting operation



Begin with step 2 if you are already in the settings menu.



1. Activate the "Menu" button.



2. Activate the menu dialogue Operation.



3. Select the desired operating mode (e.g. Therapy).



4. Confirm selection with

☞ 0K.

The device changes to therapy mode.

The **Demo** operation is only for demonstration purposes. Treatment is not possible. the **Service** operation is only accessible to the customer service. Treatment is not possible.



In this screen, you can also switch off the device, Switch off.





 Settings
 Logbook
 RF foot-switch

 Volume
 Therapy
 PIN 0 0 0 0

 Options
 Demo

 Operation
 Service

 Language
 1 2 3 4 5

 Time
 6 7 8 9 0

 PIN
 OK





if you are in the menu – settings – operation – demo and would like to change to therapy operation, the device will start the PIN request once again.

Activate Therapy.

The PIN input is blended in Switch off is blended out.

A number bar is blended in.

Enter the PIN, to do this activate the corresponding figures on the number bar.

The cursor is displayed by a green mark. On the entry of numbers, the green 0 is replaced by *.

3. Confirm the input with

SOK.



Setting the language



English


Français

Español

Begin with step 2 if you are already in the settings menu.



1. Activate the "Menu" button.





3. Select the desired language (Sz.B. English).

A Polish version is also available from SW – Version 1.3.1.



4. Confirm selection with \mathcal{S} OK.



Setting the date and time



Begin with step 2 if you are already in the settings menu.



2. Activate the menu dialogue

1. Activate the "Menu" button.



The time display is active and displays green. The cursor is displayed by a green marking of the time/date figures. A number bar is blended in.



3. Enter the time and date, to do this activate the corresponding

figures on the number bar.

If you only want to change a certain number,

activate the number (the first digit changes colour from white to green) and activate the

corresponding figures on the number bar.



4. Confirm selection with

SOK.



Changing the PIN







Begin with step 2 if you are already in the settings menu.

- 1. Activate the "Menu" button.
- **2.** Activate the menu dialogue *PIN*.
- 3. Activate the function

Change PIN.

The text *Change PIN* is replaced by the text *Enter old PIN*.

The cursor is displayed by a green mark. A number bar is blended in.

 Enter the old PIN (e.g. 1234) to do this, activate the corresponding figures on the number bar.

On the entry of numbers, the green 0 is replaced by *.

If the identification is successful, the text *Enter old PIN* is replaced by the text *Enter new PIN*.







5. Enter the old PIN (e.g. 5432) to do this, activate the corresponding figures on the number bar.

The cursor is displayed by a green marking of the PIN code.

6. Confirm the input with

☞ OK.




Setting the frequency of the PIN - request



Begin with step 2 if you are already in the settings menu.



1. Activate the "Menu" button.





3. Activate the function $\stackrel{\bigcirc}{\sim}$ *PIN request.*

 Select the desired PIN request, to do this, activate the corresponding selection (e.g. adily).



The device emits a warning in case of selection of the PIN request *Daily* or *Off.* Read carefully and follow the warning text on the screen of the device.



- 5. Confirm the input with
 - SOK.

The device changes to therapy mode.

Reading the logbook











Begin with step 2 if you are already in the settings menu.

- 1. Activate the "Menu" button.
- 2. Activate the function

S Logbook

The Logbook display changes colour to green and moves to the left. The menu dialogue Operation is active and displays green. You can read the information concerning operating hours and laser hours on the screen.

3. Activate the menu dialogue

Therapy.

You can read the therapy data on the screen. You can "scroll" between entries by activating the arrows.

4. Activate \mathcal{S} OK to return to therapy mode.



Triggering the therapy beam

⚠ WARNING

Warn all persons in the laser area before triggering the laser beam.

All persons present in the laser area must wear safety goggles.

⚠ WARNING

The laser button may not be pressed until the fibre optic conductor or application handpiece is directed to the field of operation.

⚠ WARNING

Contact cutting involving high pressure on the fibre can lead to breakage of the fibre tip.

Avoid high pressure of the fibre to tissue.

⚠ WARNING

The fibre optic conductor may not be unplugged during the therapy.

⚠ WARNING

Monitor the tissue effect continually and visually during the therapy.

⚠ WARNING

It is essential to follow the instructions for use of the fibre optic conductor.



if an error occurs during the therapy, the device will show an ERROR number and an explanatory text. In this case carry out a new start. If the error persists, please notify the KaVo Service and tell the service employee the error message, the software version and the serial number of the device.



Triggering the therapy beam with the wired foot switch

The laser changes to "Standby" as soon as the fibre optic conductor and the foot switch are connected and there is no disturbance.

- 1. Select the indication.
- 2. Check or set the therapy parameters.
- 3. Activate the laser bar or the switch-on button of the device.



The laser bar "Standby" is replaced by the laser bar "Ready". On transition, the background colour changes from yellow to green within 2 s from the centre to the outside. Hence the foot switch is released for laser activation.



 Activate the foot switch when the laser bar "Ready" has completely changed the background colour to green.

During laser radiation, the laser bar with the display LASER is displayed flashing white/red and a laser warning signal is audible. The date and time are not displayed in laser operation.

⚠ WARNING

Risk of injury!

Activate the laser bar or the switch-on button of the device and switch the device to "Standby" - mode if you want to interrupt the laser operation or treatment.




Triggering the therapy beam with the wireless foot switch

- 1. Connect the fibre optic conductor.
- Log the wireless foot switch into the device (see Assemble wireless foot switch and log into the device, page 55).

The laser changes to "Standby" as soon as the fibre optic conductor and the wireless foot switch are connected and there is no disturbance.


- 3. Select the indication.
- 4. Check or set the therapy parameters.



5. Activate the left button on the wireless foot switch once.



The laser bar "Standby" is replaced by the laser bar "Ready". On transition, the background colour changes from yellow to green within 2 s from the centre to the outside. Hence the foot switch is released for laser activation.



6. Activate the middle button on the wireless foot switch when the laser bar "Ready" has completely changed the background colour to green on the touchscreen.



During laser radiation, the laser bar with the display LASER is displayed flashing white/red and a laser warning signal is audible. The date, time and wireless foot switch symbol are not displayed in laser operation.

⚠ WARNING

Risk of injury!

Activate the left button on the wireless foot switch and switch the device to "Standby" - mode if you want to interrupt the laser operation or treatment.




Switching off the device



Every laser device in current disuse must be safeguarded against unauthorised use.



Begin with step 2 if you are already in the settings menu.



1. Activate the "Menu" button.



2. Activate the menu dialogue Operation.



3. Select Switch off.



4. Confirm selection with OK.

The device switches itself off.

⚠ CAUTION

Pull the power plug to disconnect the device from the power supply.



Transportation of the device

⚠ CAUTION

Intense vibration cab cause defects in the device that mandates a safety check.

Do not expose the device to excessive vibration during transportation.

⚠ CAUTION

Damage to the device during transportation is possible.

Take measures to prevent damage to the screen and cladding.

- Switch the device off.
- **2.** Disconnect the power plug from the power supply.
- 3. Disconnect the foot switch cable from the device.
- **4.** Disconnect the door contact plug if existing.
- 5. Roll the fibre optic conductor onto its holder.
- 6. Carry the device as shown below.

The manufacturer recommends transporting the wireless foot switch by the handle (see also the note on page 38).









⚠ CAUTION

To transport the device and accessories, use the case (see 0 to Figure 9, from page 33).



Maintenance, cleaning, safety check



Maintenance carried out by unauthorised persons could lead to life-threatening personal injury and/or severe material damage to the device.

Maintenance may only be performed by authorised persons. Authorised persons are exclusively persons who have been trained by KaVo or by a company authorised by KaVo. This group of persons are provided with the documents necessary for the maintenance and servicing. Authorised persons for device maintenance are e.g. KaVo Service staff.



The maintenance of accessories is described in the instructions for use of the accessories. The instructions for use of accessories must be observed.

Maintenance by the operator

⚠ CAUTION

Prescribed routine maintenance work may only be performed by trained and authorised practice and clinical personnel.

Table 6. Prescribed routine maintenance work

Maintenance	Interval
Please check and clean the outer surfaces	before/after use
Check the accessories	before/after use
Check the foot switch and the connection to the device	before/after use
Check the laser warning light on the OP door	before/after use
Check the electric cable (power cable and accessory cables)	before/after use
Check the functioning of the LASER-STOP button	before use
Check the battery charging status of the device and wireless foot switch (MASTER lase)	before use



Cleaning



The device must be disconnected from the power supply prior to cleaning.

Device parts, in particular electrical parts may not be exposed to water spray.

The device may not be disinfected with gas and/or spray disinfectants. Information from manufacturers of disinfectants and legal regulations with regard to disinfection and protection against explosion must be observed.

↑ CAUTION

The device may only be wiped down with a moist cloth on the outside with conventional clinical cleaning agents.

The control panel (touchscreen) must be cleaned especially carefully and gently using a dry, soft cloth. If the control panel (touchscreen) is really dirty, it can be cleaned with a slightly moistened cloth.

- 1. Switch the device off.
- Remove the fibre optic conductor from the device if 2. necessary.
- 3. Dismantle the transmission fibre optic conductor and unscrew the plug of the transmission fibre optic conductor (see Figure 15, page 107).
- Clean the device and the dismantled parts. It is essential 4. to observe the warning instructions stated above.
- 5. Carry out the wipe disinfection of the transmission fibre optic conductor.




The cleaning, disinfection (or sterilisation) of accessories is described in the instructions for use of the accessories.

The instructions for use of accessories must be observed.



Figure 15. Dismantling the transmission fibre optic conductor winding





6. Assemble the device in reverse sequence.



Safety checks

Safety checks are understood as maintenance conducted at regular intervals.



Safety checks carried out by unauthorised persons could lead to life-threatening personal injury and/or severe material damage to the device.

Safety checks may only be performed by authorised persons. Authorised persons are exclusively persons who have been trained by Kaltenbach & Voigt GmbH tor by a company authorised by Kaltenbach & Voigt GmbH.



The comparison of the laser performance at the fibre end with the value set on the touchscreen is a mandatory constituent of the safety check.



The laser performance can be checked with a conventional, calibrated laser performance measuring device suitable for the wavelength.

The interval for the device is determined at every 500 operating hours or every 12 months, at the latest. A note appears in the device display if the interval is exceeded.

After every change of location outside the clinic or dental practice - in case of appropriate transportation without detectable transport damage - a function check is sufficient (National German Guidelines / BGV A3, MPBetreibV).

Function check comprises commissioning with device self test (see page 57), fibre optic conductor with pilot laser (see page 65) and subsequent check of laser performance at fibre end at 5 Watt setting.



Messages

Messages are indicated in the display. Some messages are accompanied by an acoustic alarm.

Table 7. Messages

Message	Meaning	Action, reference
Connect fiber to device	The transmission optic fibre conductor is not or not correctly connected	Then check whether the transmission optic fibre conductor is correctly connected.
Unlock remote interlock	The door contact is interrupted	Check whether the door contact switch is closed.
Over- temperature housing	the device has switched off because the housing temperature or therapy laser beam is too high	Wait until the device is ready for operation again (do not switch off). Please contact KaVo Service if the message is still displayed after 2/5 minutes.
Over- temperature fiber connector	The device has switched off because the temperature at the transmission fibre optic conductor or therapy laser beam is too high	Check the proximal end of the transmission fibre optic conductor (with a magnifying glass). If you do not detect any damage, wait until the device is ready for operation again (do not switch off). If you detect any damage when you check the transmission fibre optic conductor, replace the conductor. Please contact KaVo Service if the message is still displayed after 2-5 minutes.
Over- temperature laser	The device has switched off the therapy laser beam because the laser diode temperature too high	Wait until the device is ready for operation again (do not switch off). Please contact KaVo Service if the message is still displayed after 2/5 minutes.
Over- temperature heatsink	The device has switched off the therapy laser beam because the cooling unit is overloaded	Wait until the device is ready for operation again (do not switch off). Please contact KaVo Service if the message is still displayed after 2/5 minutes.
High temperature, switch-off possible	The battery temperature is too high.	The current therapy can be continued for 3 - 5 minutes.



Message	Meaning	Action, reference
Over- temperature battery pack	The battery temperature is too high.	Do not switch off the device. Wait until the device is ready for operation again. The therapy can be continued in approx. 30 minutes.
Heatsink non- functional	The cooling unit is defective. The device switches off the therapy laser beam	Switch off the device, please notify the KaVo Service.
Deviation laser power	The deviations in the performance of the laser diodes are outside the tolerance limits. The device switches off the therapy laser beam.	Continue the therapy. Please contact KaVo Service if the message is displayed repeatedly.
Pilot control not possible	The pilot diode has overtemperature or is defective	Wait until the device is ready for operation again (do not switch off). Please contact KaVo Service if the message is still displayed after 2/5 minutes.
Do not tilt RF foot-switch	The wireless foot switch is considerably tilted for longer than 1.5 seconds	If this position of the foot switch is not rectified within 5 sec., the foot switch logs off automatically. You must log the wireless foot switch into the device again to continue the therapy. Put the foot switch in a horizontal position and continue the therapy. Replace the foot switch if the message appears repeatedly.
Update date/time	Invalid date and time set	Set the correct date and time.
Exchange RTC battery	Index	Notify service, have the battery for the real/time clock changed.
Battery charging insufficient	Device battery charging is lower than 25%	Charge the battery completely. If you are administering therapy, connect the device to the power supply and charge the battery completely after the therapy.
Battery health weak	The battery has reached the end of its life	Notify service, have the battery changed.



Product service life and disposal

The service life of the instrument is designed for an application period of at least 8 years.

At the end of the service life of the laser system, contact the manufacturer or his representative to obtain instructions for disposal or recycling (WEEE 2002/96/EC: European Guideline for the Disposal of Obsolete Electrical Devices).

At the end of the service life cycle, dispose of the optic fibres (fibre optic conductors) in accordance with the local regulations for contaminated products. In this way you will minimise hazards for the environment and personnel that could be caused by contamination and residues of the fibres employed.

Disposal of accessories in accordance with the instructions for use of accessories.

China-RoHS declaration of conformity for products of the company Kaltenbach & Voigt GmbH

This product of the company Kaltenbach & Voigt GmbH is complaint with the Industrial Directive # 39 of the Ministry of Information of the People's republic of China, adopted on 28. February 2006, coming into force on 1. March 2007 ("Control of Environmental Pollution by Electronic Products").

This product bears the EFUP environmental protection mark (environmentally-friendly period of usage), in accordance with the industrial standard SJ/T11364-2006.

This mark indicates that while the product contains specific toxic or hazardous substances, it does not cause environmental pollution and can be used safely during its EFU period. This period is defined for a specific time period and the corresponding mark is attached to the product (e.g. 50 = 50 years EFUP).

Figure 16. Environmentally-friendly service periods – marking/China

