GARANTIEZERTIFIKAT

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GUARANTEE CERTIFICATE

FARO offers the final customer a **24 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser. **Repair under guarantee is considered valid only when:**

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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SYMBOLS

DANGER

The paragraphs marked with this symbol contain instructions that must be carried out carefully to prevent damage to the system, the operator, and possibly the patient.

WARNING

These instructions advise that great care must be taken to avoid situations that could damage the system.

PROHIBITION

This icon highlights what should not be done in order to prevent damaging the device.

SUGGESTIONS

With this icon, information is provided that enables the system to be used more efficiently.

Dispose of the system in accordance with regulations for the separate collection of electrical material.

Device classe II

Consult the enclosed documentation.

SAFETY REGULATIONS

- The function of the ALYA dental light is to illuminate the oral cavity of the patient during dental operations.
-  - The light must be used by qualified staff.
- The articulated arm and head joint ensure that the light beam is correctly positioned. Do not overload the arms and joints with impact on the limit switches.
- Installation of the system must only be performed by specialized staff.
- Check that the power voltage, shown on the data plate, corresponds to the mains voltage.
- Do not block the light beam for patients at risk (i.e. with eye conditions). Use correct protection and precautions.
- Information plate with LED indications fixed to the rear arm.
- The dental light must be installed on a specific control and supply device, such as a dental unit, or **with an electrical system that complies with standard IEC 364-1 and "national regulations for installation of electrical systems in areas for medical use".**
-  - The device must be installed with an omnipolar shut-off switch that disconnects it from the mains and is conformant to IEC/EN 61058. It is necessary to install a green signal light to indicate when the lamp is switched on.
-  - Do not carry out any maintenance work when the power is on; disconnect the power cable from the mains beforehand.
-  - Do not introduce items or tools into slits in the light head that could touch live points.
- Patients at risk (e.g. children with eye diseases) must not stare at the light beam. Always use appropriate protection devices and precautions. Faro suggests using the protective eyewear BLUE-BAN or BABY BLUE-BAN (for children).
-  - Not suitable for installation in environments containing flammable gas.
- Do not spray detergents-disinfectants directly on the head.
- Inadequate performance of the system does not in any way compromise the safety of the patient.
- For cleaning plastic parts of the lamp ALYA, **DO NOT USE** detergents-disinfectants containing: **AMMONIUM HYDROXIDE – SODIUM HYDROXIDE – METHYLENE CHLORIDE – METHYL ALCOHOL**. Failure to comply with this indication might cause: • risk of breakage of plastic parts • risk of structural breakdown of joints with potential fall of the articulated arm. In case of doubt, please contact CUSTOMER CARE FARO.
-  - Do not carry out maintenance or replacement work on parts other than those shown in the manual. Any operation not indicated in the manual could compromise the safety aspect of the device.
- Disinfect surfaces with hydroalcoholic disinfectants.
- Do not leave small parts unattended or within the reach of exposed persons (children); they could be potential sources of danger.
- The dentist is obliged to use the disposal protection sleeves on the handles of the lamp or ensure that versions without this possibility are sterilized.
-  - Divide materials according to their type (e.g. iron, rubber, plastic.). When scrapping and disposing of the materials, comply with the current legislation of your country, and if necessary use recognized and authorized specialist companies.
-  - The packaging of the light is suitable for protecting it from the penetration of external agents.
- The system in its original packaging can be transported or kept in storage for a period of 15 weeks if the following environmental conditions are observed:
 - **Ambient temperature from -20° C to + 70°C**
 - **Relative humidity from 10% to 90%**
 - **Atmospheric pressure from 500 to 1060 mBar**
- The system must be used in the following environmental conditions:
 - **Temperature from 10° to 40° C**
 - **Relative humidity from 30 to 75%**
 - **Atmospheric pressure from 700 to 1060 mBar**

SAFETY REGULATIONS

REQUIREMENTS FOR ELECTROMAGNETIC COMPATIBILITY

This medical device requires particular precautionary measures to ensure electromagnetic compatibility, and must be installed and used in compliance with the information provided in the accompanying documentation.

Manufacturer's guide and declaration – Electromagnetic emissions		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Emission tests	Conformity	Guide – Electromagnetic environment
RF Emission CISPR15	Yes	The ALYA lamp uses RF energy only for internal functioning. Therefore, its RF emissions are very low and probably do not cause any interference with other nearby electronic equipment.
RF Emission CISPR15	Yes	The ALYA lamp is suitable for use in all buildings, including domestic buildings and those connected directly to the low voltage public electricity mains that power domestic buildings.
Harmonic emission	Class C	
RF Emission CISPR11 / EN 55011	Yes	The ALYA lamp is not suitable for interconnections with other devices (ceiling model).

ELECTROMAGNETIC IMMUNITY

Manufacturer's guide and declaration – Electromagnetic immunity		
The ALYA lamp is made to function in the electromagnetic environment specified below. The client or user should make sure that it is effectively used in this environment.		
Immunity test	Conformity	Guide – Electromagnetic environment
Electrostatic discharge (ESD) IEC/EN61000-4-2	$\pm 6\text{kV}$ contact $\pm 8\text{kV}$ air	Flooring must be in wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC/EN61000-4-4	$\pm 2\text{kV}$ power supply $\pm 1\text{kV}$ for input/output lines	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Surge IEC/EN61000-4-5	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment.
Voltage dips, short interruption and voltage variation IEC/EN61000-4-11	$< 5\%$ Ut for 0,5 cycle 40% Ut for 05 cycle 70% Ut for 25 cycle $< 5\%$ Ut for 5 sec.	The quality of network voltage should be equivalent to that of a typical commercial or hospital environment. If the user of the ALYA lamp requires continual use, also when electrical power is not present, it is advisable to have a UPS (uninterrupted power supply).
Power frequency magnetic field IEC/EN61000-4-8	3A/m	Level of the magnetic field at the typical power main frequency of a commercial or hospital environment.
Conducted immunity IEC/EN61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting equipment)	Portable and mobile RF communication devices should not be used near any part of the dental unit, including cables, except for when they comply with the recommended separation distances calculated by the equation applicable for the transmitter frequency. Recommended separation distances: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ from 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ from 800 MHz to 2,5 GHz Where P is the maximum nominal power emitted from the transmitter in Watts (W), according to the transmitter manufacturer, and d is the recommended separation distance in metres (m). The intensity of the fixed RF transmitter field, as determined in an electromagnetic study of the site a, could be less than the level of conformity of each frequency interval. Interference could be experienced in proximity to devices bearing the following symbol: 
Conducted immunity IEC/EN61000-4-6	3Vrms 80MHz to 2.5GHz (for non life-supporting equipment)	

Note: Ut is the value of the voltage.

Note 1: At 80 MHz and 800 Mhz the higher frequency interval is applied.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.

a) The ISM bands (industrial, scientific and medical) between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz.

b) The levels of conformity of the ISM bands between 150 kHz and 80 MHz and in the bands 80 MHz to 2.5 GHz are indicative. The probability of a portable transmitter causing interference if accidentally brought in the range of a patient decreases depending on environment conditions. This is why an additional factor of 10/3 has been incorporated into the formula used to calculate the separation distance from transmitters.

c) Field intensities for fixed transmitters like base stations for radio-telephones (cellular and cordless phones) and terrestrial transmitters, amateur radio stations, AM and FM radio transmitters and TV transmitters cannot be predicted theoretically and precisely. To establish an electromagnetic environment caused by fixed RF transmitters, an electromagnetic investigation of the site should be done. If the intensity of the field measured in the area of use of the dental unit is higher than the applicable level of conformity described above, it is necessary to monitor the lamp to check that it functions properly. If abnormal performance is detected, additional measures could prove necessary, such as a different orientation or position of the lamp.

d) The field intensity in an interval of frequencies from 150 kHz to 80 MHz should be less than 3 V/m.

**Recommended separation distances between portable and mobile
radio-communication devices and the dental unit**

The ALYA lamp is designed to function in an electromagnetic environment in which disturbances from RF emissions are kept under control. The client or operator of the unit can contribute to preventing electromagnetic interference by assuring a minimum distance between mobile and portable RF communication devices (transmitters) and the dental unit, as recommended below, in relation to the maximum output power of the radio-communication devices.

Maximum nominal output power of the transmitter W	Separation distance at the transmitter frequency m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal power emission not listed above, the recommended separation distance d in metres (m) can be calculated using the equation applied to the transmitter frequency, where P is the maximum nominal power emission of the transmitter in Watts (W) according to the manufacturer.

Note:

At 80 MHz and 800 MHz, the highest interval is applied.

These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and people.


FEATURES

Versions

The "ALYA" dental light is supplied in the following versions:

- **S/TS Light (with/without transformer) with joystick**
 - **S/TS Light (with/without transformer) with switch at the table keyboard**
 - **S/TS Light (with/without transformer) with proximity sensor**
 - **S/TS Light (with/without transformer) with ceiling mounted version with joystick or proximity sensor**
- **Light source** formed by two LEDs whose light is mirrored on two dishes, passing through 2 secondary lenses.
 - **Reflecting surfaces** with dishes that allow to obtain a spot of regular uniform light at all lighting levels and to uniformly distribute the light in the operating field without the operator creating shadows or dark areas.
 - **Light intensity adjustment** using a mechanical switch or with a sensor.
 - **Infra-red sensor** enables light to be switched on or off without direct contact, thereby eliminating the possibility of cross-infection.
 - **Easy maintenance** thanks to the application of new technology that takes into consideration various requirement in terms of safety, ergonomics and hygiene.
 - **Removable handles** for sterilization purposes.

DESCRIPTION OF PARTS

- 1 - Joystick
- 2 - Lens
- 3 - Parabola
- 4 - Handle
- 5 - Proximity sensor
- 6 - Fuses

INSTALLATION AND CONNECTIONS**“ALYA” Dental Light, version S/TS**

Check that the package contains the following components:

- Dental light (in the version requested)
- Lever envelope + tool
- Instruction manual

The system must be cleaned before use.

(see the paragraph Cleaning the device)

The system must be installed by specialized technicians.

The power must always be off when the system is installed.

Simultaneous use of the lamp with an electric scalpel can cause the scalpel to malfunction. To reduce risks, it is advisable to use the version controlled with the joystick instead of the proximity sensor. Always verify compatibility before use.

Lamp assembly, combined version

- Install the lamp by inserting the terminal lamp pin into the specific hole on the combined unit.


	SAFE WORKING LOAD	MINIMUM BREAKING LOAD
Arm length 855 mm	3.12 Kg	26 Kg.
Arm length 550 mm	2.64 Kg	35 Kg

1) POWER CABLE:

- **LIGHT WITH TRANSFORMER;** must be powered directly from the mains, which should correspond to the specifications indicated on the data plate or in the technical specifications of the manual.
- Ensure that the line switch conforms to the surface distances in air as specified in International Publication IEC 328.
- **LIGHT WITHOUT TRANSFORMER;** must be powered with low-voltage alternating current (24V AC) using a safety transformer that conforms to EN 60601-1.

2) REMOTE CONTROL CABLE (where available):

- Connect the cable with two push-buttons (not supplied), whose contact is normally open, according to the following scheme:

INSTRUCTIONS FOR USE

SWITCHES SYMBOL

ON-OFF

“ALYA” LIGHT WITH JOYSTICK

(see the description of parts)

Turning on and off/ Adjustments:

- To turn the lamp on and off, press and release the joystick lever to the left or right side.
- Adjustment:
 - a) to reduce light intensity, keep the joystick lever pressed to the left side (rear view of lamp) until it reaches the desired intensity. When it reaches the minimum intensity, the machine will beep once.
 - b) to raise the light intensity, keep the joystick lever pressed to the right side (rear view of the lamp) until it reaches the desired intensity. When it reaches the maximum intensity, the machine will beep twice.

The control joystick must be handled delicately to prevent breakage.

Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.

“ALYA” LIGHT WITH PROXIMITY SENSOR

(see the description of parts)

Turning on and off:

- To switch on and off, move near to the sensor up to a maximum distance of 3 cm
- For **light intensity** regulation, place the hand near the sensor until desired intensity is reached, from the maximum to the minimum level and from the minimum to the maximum level. On reaching maximum intensity, an acoustic signal will be heard (2 beeps); there will be 1 beep for minimum intensity.

Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.

“ALYA” DENTAL LIGHT SYSTEM WITH REMOTE CONTROL

(see the description of parts)

Turning on and off/ Adjustments:

- To switch on/switch off the lamp, press the push-button “A”.
 - Adjustment:
 - a) To reduce the light intensity keep the push-button “A” pressed, until the desired intensity is obtained.
When the minimum light intensity is obtained, you will hear an acoustic signal (1 beep).
 - b) To increase the light intensity keep the push-button “A” pressed, until the desired intensity is obtained.
When the minimum light intensity is obtained, you will hear an acoustic signal (2 beeps).
 - c) To immediately obtain the minimum light intensity, press the push-buttons “B”.
- When the minimum light intensity is obtained, you will hear an acoustic signal (1 beep). Press the button again to restore previously selected light intensity.


Each time the lamp is turned on, the light intensity will be at the level memorised when it was turned off the time before.

ALYA SWITCH LEVER ASSEMBLY

- Insert lever "A" on switch pin.
- Place lever "A" hole in correspondance of "B" pin axe.
- Screw dowel "C" with the Allen wrench provided.

 The control joystick must be handled delicately to prevent breakage.

**MAINTENANCE****Replacing the fuses "6" in the version with transformer (see the description of parts)**

The light with transformer has two fuses of the same value.

They are replaced in the following way:

- Make sure that the power is off.
- Unscrew the plugs "6" on the transformer holder.
- Remove the fuses. After verifying the problem, replace them.

 It is important that the new fuses have the specifications given on the data plate and in the technical specifications.

CLEANING THE PARABOLAS "3" (see the description of parts)

For cleaning, Faro recommends using the detergent "Faro Perfex"; as an alternative, clean with a soft cloth and isopropyl alcohol. Do not use detergents containing surfactants or water repellent agents, whose build up can leave halos.

 Hydroalcoholic disinfectants with 70% isopropyl or ethyl alcohol are suitable for disinfection.

For cleaning plastic parts of the lamp ALYA, DO NOT USE detergents-disinfectants containing **AMMONIUM HYDROXIDE - SODIUM HYDROXIDE - METHYLENE CHLORIDE -METHYL ALCOHOL**.

Faro has tested and suggests the following disinfectants: • Eco Jet-1 (Cattani Group) • Sporekin Plus DS (Ims srl) • Zefiro Quick (Molteni Dental). After disinfection, Faro recommends using the product "Faro Perfex" for removing residue and halos.

 PLEASE NOTE: other products might damage the reflectors (dishes). In case of doubt, please contact CUSTOMER CARE FARO.

 Slight halos do not impair the quality of light.

"A"

**STERILIZING THE HANDLES**

To remove the handle, unscrew button "A" and slip it off. To insert it, push firmly and screw on "A".

 Handles are not provided sterile and must, therefore, be sterilised before use.

Sterilize the device with standard cycles (SHPGRP) at 121°/134° C for a total of 200 sterilisation cycles.

OTHER PARTS OF THE LIGHT (head-articulated arm)
Clean with a soft cloth.

 It is strictly forbidden to use abrasive substances, detergents with trichloroethylene, benzene, oil of turpentine or similar on any parts of the light.

 Periodically perform the electrical safety check recommended for devices installed in medical premises.

PERIODIC CHECKS

- Check the absence of space between joint components of the arms (**yearly**)
- Check readability of data on the information plate (**yearly**)
- Power safety check: (**every two years**)
 1. Rigidity
 2. Dispersion
- Light checks (**every five years or 10,000 hours of function**)
 1. Maximum lighting intensity >35000 lux
 2. CRI decay: <20% for data on the information plate.
 3. Expected value of Blue light on spectrum emitted, measured in w/m²: <100.

INSTALLING THE CEILING LIGHT FIXTURE

Before installing the lamp on a structure (e.g. ceiling), calculate the overall load of the system to be installed by adding the weight of the system to the maximum load that can be applied to the lamp. Then, ascertain that the chosen structure and fixing elements are suitable for bearing the resulting load multiplied by a safety factor K8.

Also ascertain that fixing elements are suitable for ceiling materials (e.g. concrete, stone, solid brick, hollow brick).

For ceiling-mounted FARO light fixtures designed to hold 2 lamps, calculate the load as specified below:

- Alya lamp without transformer and complete with arm and hinge joint:	6 kg
- Column for ceiling-mounting, length 1700 mm:	8 kg
- Ceiling lamp with 2 transformers for application of 2 lamps:	8 kg
- Load applicable to the lamp:	35 kg
- Total load for 2 lamps: $6 \times 2 + 8 + 8 + 35 \times 2 =$	98 kg

Load with safety coefficient 8: 98 kg X 8 = 784 kg (7.7 kN)

SUPPLIED

N° 5 - Expansion bolts (ideal for concrete and natural stone)

ACOUSTIC SIGNALS

OTP*	= Beep 30 seconds
MAX	= 1 beep - 2 beep with proximity
MIN	= 1 beep
1 BEEP	= On the commands
1 BEEP	= At start-up
OpL**	= Beep 30 seconds

* OTP: Over-temperature protection

** OpL: LED charger unplugged

TROUBLESHOOTING

The light does not switch on

- Check that the power supply is on.
- Check the state of the fuses.
- If none of these causes is verified, consult technical assistance.

The light intensity is considerably reduced

- Clean the parabolas.
If the light intensity does not return to initial values, consult technical assistance.

The handles do not hook on, or they are difficult to hook onto the machine.


- Carry out a sterilisation cycle at 121°/134° with the specific core former PN 310136.

Patches have appeared on the reflectors (dishes)

- Clean the surfaces with the specific product "Faro Reflex".
- Clean the surfaces with isopropyl alcohol.
- An unsuitable product for cleaning or disinfection has been used that has damaged the surfaces.
- Contact customer care FARO.

TECHNICAL SPECIFICATIONS

Power voltage (without transformer)	: 17÷24V ac ± 10% - 50/60 Hz 22÷35V dc
Power voltage (with transformer)	: 230 V 50/60 Hz
Absorbed power	: 26 W
Fuses (version with transformer)	: 2 x T250mA 250V
Protection against electrical hazards	: Device class II

Label compliant with EN 60825-1:94+A1:02+A2:01

OPTICAL SPECIFICATIONS IN COMPLIANCE WITH ISO 9680

Luminous spot dimensions	: 180 mm x 90 mm
Lux	: 3000*-50000* lux @700mm
Colour temperature	: 5000 K
CRI (color rendering index)	: >95*

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