

KaVo Everest® portal.



Increase your efficiency!

The Everest Portal is a network that ensures an exchange of 3D-construction data and therefore acts as an interactive "exchange-platform" for several scan and milling partners.

What are the advantages for the portal?

- Increase the efficiency of Everest laboratories by optimizing capacity
- Excellent marketing-tool – extension of your customer basis!
- Interactive exchange platform for you and numerous scanning partners
- Extended network capability and use of various features of Everest-systems – connects Everest users worldwide
- Design control remains in the laboratory and not in the milling center

Which conditions are necessary for the Everest portal?

- Everest Scan Pro, Everest Scan II, 3shape Scanner
- high speed internet
- latest software version
- Windows XP

How does the Portal work?

- The milling laboratory sends a manufacturing request into the KaVo Everest portal
- The Everest scan laboratory is now able to see all registered milling partners within the portal which have free capacities available. Hence the laboratory can choose a suitable partner.
- After accepting the offer construction data will be transferred securely and without complications to the milling partner by the Everest portal (a KaVo server).
- The milling laboratory manufactures the order and sends it back to the scan laboratory.
- Writing and payment of the invoice take place directly between the milling and the scanning laboratory.
- There is also the production center BEGO Medical in the portal, available for the manufacturing of laser-sintered NEM work.



You will find the portal software on www.kavo-everest.com for downloading.




KaVo. Dental Excellence.

KaVo Everest portal

Through the Everest portal you can offer free capacities and therefore your Everest System will be utilised ideally and economically. Also construction data of the scan laboratories will be passed quickly and easily by the portal.

The KaVo Everest portal as a worldwide exchange-platform for scan and milling laboratories will increase your profitability. Therefore you will use the variety of the Everest systems ideally!



Interested?

Visit us on our homepage www.kavo-everest.com and get into the variety of the Everest-community by downloading the portal-software for free.

Our hotline will be glad to give you advice! Phone +49 73 51 56-22 00 or -21 00



KaVo. Dental Excellence.