Sella II **SEALER**



AND INNOVATION

Technical Features

AC 110/220 V 50/60 Hz **Power input Power otput**

100 W 12 mm.

30 cm.

Pre-heating time **Heat preservation** Yes

Dimensions (W*L*H) 46*36*17 cm **Gross weight** 7.5 Kg.

Sella II

Seal width

Seal length

is intended for sealing envelopes for steam sterilization, one of sterilization auxiliary series products, which is attractive, stable, easy operation and low maintenance rate. Terrific transformer will automatically adjust to the best sealing temperature even in unstable voltage. Nice aluminum alloy body, easy for clean, and come with integrated cutter with 100,000 times cutting guarantee

- Sealing length: 30 cm Sealing width: 12 mm
- According to European standard sealing
- Bi-directional operational handle design
- PC board with automatic temperature control, available for continuous using
- Beep-alert (for heating Status)
- 100W energy saving design
- Comes with integrated cutter with 100,000 times
- cutting guarantee and roll support



via Faro, 15 - 20876 Ornago (MB) - Italy Tel. +39 039.68781 - Fax +39 039.6010540 www.faro.it - export@faro.it