Start safe with affordable KaVo quality.

KaVo smart

series

KaVo SMARTtorque

- Price first. KaVo Standards at a modest price and with the flexibility of different couplings
- Safety first. KaVo chucking system with up to 30 N chuck retention force and 4-port spray
- My first KaVo. Enjoy KaVo quality by the use of hard ceramic ball bearings in standard and mini head versions



KaVo SMARTtorque LUX S619 L

Mini LUX S615 L



KaVo SMARTtorque S619 L/C / S615 L Mini

- Light (L) and non-light (C) versions for MULTIflex couplings
- Head versions: standard (h: 14.5 mm) or mini (h: 12.3 mm)
- Reliable 16 Watt power (Mini: 13 Watt), for efficient and safe preparation
- High safety, with chuck retention force up to 30 N
- 4-port spray (Mini: single spray), for optimum cooling effect
- Enhanced, anti-reflux effect
- New long-life, ceramic bearings

KaVo SMARTtorque S609 C / S605 C Mini

- 4-hole standard connection (Midwest)
- · Head versions: standard (h: 14.5 mm) or mini (h: 12.3 mm)
- Reliable 16 Watt power (Mini: 13 Watt), for efficient and safe preparation
- High safety, with chuck retention force up to 30 N
- 4-port spray (Mini: single spray), for optimum cooling effect
- Enhanced, anti-reflux effect
- New long-life, ceramic bearings

KaVo SMARTtorque S608 C

- 2-hole standard connection (Borden)
- Head version: standard (h: 14.5 mm)
- · Reliable 16 Watt power, for efficient and safe preparation
- High safety, with chuck retention force up to 30 N
- 4-port spray, for optimum cooling effect
- Enhanced anti-reflux effect
- New long-life, ceramic bearings



Particular SMARTtorque Features:

- Low price for entry into the KaVo world of dental excellence
- 12 months warranty
- High torque and perfect cooling spray for rapid, efficient preparation
- KaVo safety standards, such as a chuck retention force up to 30N
- KaVo Push-button chuck system, for rapid bur changing
- 6 different models / price versions

Made in Germany:

- Proven quality, that has made KaVo into the worldwide market leader
- 100 years of experience goes into each product
- High quality preparation, such as one-piece handpiece casings

KaVo smart

Item	Description	Mat. No.
S619 L	SMARTtorque (Light) with MULTIflex connection	1.008.1641
S619 C	SMARTtorque (non-light version) with MULTIflex connection	1.008.1642
S615 L Mini	SMARTtorque Mini (Light) with MULTIflex connection	1.008.1643
S609 C	SMARTtorque with Midwest-connection	1.008.1644
S605 C Mini	SMARTtorque Mini with Midwest-connection	1.008.1645
S608 C	SMARTtorque with Borden-connection	3.001.6560
MULTIflex 460 LED	MULTIflex LUX coupling 460 LED	1.007.3201
MULTIflex 460 LE	MULTIflex LUX coupling 460 LE	1.001.7599
MULTIflex 460 E	MULTIflex coupling 460 E (non-light version)	1.001.7600







KaVo. Dental Excellence.