



Price first. Safety first. My first KaVo.



**KaVo SMARTtorque –  
The new turbine from KaVo**

Be smart and play safe with KaVo SMARTtorque turbines. For KaVo safety features at an extremely affordable price. Just another instrument from the KaVo SMART Series – the perfect chance to step into the KaVo world of quality.



KaVo. Dental Excellence.

KaVo SMARTtorque.

**Flexible, individual, low priced**

1.008.1641	SMARTtorque LUX S619 L
1.008.1642	SMARTtorque S619 C
1.008.1643	SMARTtorque mini LUX S615 L
1.008.1644	SMARTtorque S609 C
1.008.1645	SMARTtorque mini S605 C
3.001.6560	SMARTtorque S608 C
1.007.3201	MULTiflex LUX coupling 460 LED
1.001.7599	MULTiflex LUX coupling 460 LE
1.001.7600	MULTiflex coupling 460 E



**KaVo smart  
series**

Start safe with affordable KaVo quality.

» Entry level price  
for the world of KaVo

» Solid  
basic equipment

» Made in  
Germany

**KaVo expert  
series**

Rely on advanced KaVo performance.

» Best price / longevity  
ratio

» Proven  
configuration solutions

» Highest  
efficiency

**KaVo master  
series**

Enjoy maximum power with innovative  
KaVo comfort.

» Highest  
treatment quality

» Maximum treatment  
comfort

» Increased patient  
satisfaction



KaVo. Dental Excellence.

KaVo Dental GmbH · D-88400 Biberach/Riß · Telefon +49 7351 56-0 · Fax +49 7351 56-1488 · www.kavo.com

KaVo SMARTtorque.

Price first.  
Safety first.  
My first KaVo.

**KaVo smart  
series**

Start safe with affordable KaVo quality.



**NEW**

KaVo SMARTtorque  
Made in Germany



KaVo. Dental Excellence.

## Price first.

### KaVo SMARTtorque – KaVo Standards at a modest price

- Solid basic equipment, KaVo quality at a low price
- Six different models / price versions:
  - 3x SMARTtorque with MULTIflex connection:**
    - > 2x SMARTtorque with standard head, with and without light (H: 14.5 mm, head diameter: 12.5 mm)
    - > 1x SMARTtorque with mini head with light (H: 12.3 mm; head diameter: 10.8 mm)
  - 2x SMARTtorque with 4-hole standard connection (Midwest):**
    - > 1x SMARTtorque with standard head
    - > 1x SMARTtorque with mini head
  - 1x SMARTtorque with 2-hole standard connection (Borden):**
    - > SMARTtorque with standard head without light



## Safety first.

### KaVo SMARTtorque – Safe for you and your patient

- Securely held bur or instrument, thanks to the KaVo chucking system, with up to 30N chuck retention force – no loosening or slipping of the rotating instrument
- Optimum protection of the tooth from thermal damage, by 4-port spray
- Safe, unimpeded access to all quadrants in the buccal cavity, thanks to ideal head sizes:
  - > Standard version: H: 14.5 mm, head diameter: 12.5 mm
  - > Mini version: H: 12.3 mm; head diameter: 10.8 mm
- Simple operation with the KaVo, solid glass-rod light-conductor, with 25,000 LUX at the bur tip.
  - > Glare-free, bur-concentric illumination
  - > Reliable identification of preparation margins
- Secure handling and simple cleaning thanks to an ergonomic form that is optimised for hygiene



\*optional



## My first KaVo.

### KaVo SMARTtorque – Welcome to the KaVo world of quality

- Efficient and fast preparation thanks to 16 Watt output
- Compliance with exacting KaVo standards, including for example KaVo turbine cartridges, with long-lasting, ceramic bearings
- Increase in both hygiene and ceramic bearing longevity, owing to the cartridge, which incorporates a proven, automatic, contamination-minimising, anti-reflux design
- Precise and quiet running due to the highest standards of precision in production and proven KaVo technology
- Efficient operation with the KaVo push-button chucking system, for rapid bur changing
- **KaVo – Made in Germany** – stands for over 100 years of experience in manufacturing and selling over 5,000,000 turbines, straight & contra-angled handpieces. Above all, it stands for precision, reliability, longevity and the highest standards of quality

