




Cattani air care AC310

CATTANI
air technology

OIL-LESS COMPRESSOR CAD-CAM VERSION 8-10 BAR

Il compressore a secco Cattani (privo di lubrificazione ad olio) completo di impianto di essiccazione e di una serie di filtri, è in grado di fornire aria compressa secca e di elevata qualità e purezza. Con il filtro sterilizzante (opzionale) l'aria viene inoltre trattata con efficienza del 99,9999+% su particelle da 0,01 µm. Abbiamo fatto della qualità ed affidabilità la nostra bandiera, per questo siamo in grado di dare 3 anni di garanzia.

The CATTANI oil-less compressor (without oil lubrication) complete with drying system and a set of filters, is capable of supplying dry compressed air of extremely high quality and purity. With the sterilising filter air (optional) is treated with an efficiency of 99.9999+% on 0.01µm particles. We have chosen quality and reliability as our emblems: this is why we can grant you a 3-year warranty.



Tipo	Dati elettrici	Codice
AC310	(1~) 230-50 hz- 10,2 A - 1,5 kw	013310
AC310	(3~) 400-50 hz- 3,7 A - 1,5 kw	013313
AC310	(1~) 220-60 hz- 12,2 A - 1,75 kw	013305

Aria prodotta NI/min at 8 bar	165/50 Hz	205/60 Hz
Capacità serbatoio litri	45	
Peso kg	64	
Livello pressione sonora db (A)	73,6	

Type	Electrical data	Code
AC310	(1~) 230-50 hz- 10,2 A - 1,5 kw	013310
AC310	(3~) 400-50 hz- 3,7 A - 1,5 kw	013313
AC310	(1~) 220-60 hz- 12,2 A - 1,75 kw	013305

Delivered air NI/min at 8 bar	165/50 Hz	205/60 Hz
Tank capacity liters	45	
Weight kg	64	
Sound pressure level db (A)	73,6	

CATTANI S.p.A.

Via Natta, 6/A - 43122 Parma - Italy - Tel. +39 0521 607604
 SALE DEPT. FAX: +39 0521 607628 - PURCHASING DEPT. FAX: +39 0521 607855
 ACCOUNTING DEPT. FAX : +39 0521 399966
<http://www.cattani.it> Email: info@cattani.it PEC: cattani@peclife.it

Codice Fiscale e Partita I.V.A. 01720020344 - E.E.C. VAT IT 01720020344 -
 Capitale Sociale € 1.549.800,00 I.V. - R.E.A. 173616
 Registro Imprese Parma n. 01720020344

AZIENDA CON SISTEMA DI GESTIONE
 PER LA QUALITÀ CERTIFICATO DA DNV
 = UNI EN ISO 9001:2008 =